**HERPES SIMPLEX VIRUS**  
**1+ 2 IgG & 2 IgG**  
**Rapid Test**

**HERPES SIMPLEX VIRUS 1+2 IgG & 2 IgG**  
**ImmunoFLOW™ RAPID TEST**

ImmunoFLOW™ HSV Test detects IgG to HSV and HSV gG-2 (HSV type 2). The identification aids in the diagnosis of diseases caused by herpes simplex viruses and provides epidemiological information

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**Packaging: 25 test**

## HERPES SIMPLEX VIRUS 1+2 IgG & 2 IgG Rapid Test

### ImmunoFLOW™ HSV Test detects IgG to HSV and HSV gG-2 (HSV type 2)

There are two herpes antigenic types. (1,2) "Definitive diagnosis of genital herpes infections is fundamental to the management of patients and the development of strategies to prevent transmission to partners and neonates".(3) Such diagnosis has proven inaccurate when based solely on clinical history and impression.(4) Instead, virus, antigen or nucleic acid detection and classification is used for patients presenting with lesions or type-specific serological tests may be used when lesions are absent.

For type specific serology, either western blot (5,6,7) or assays or type specific protein (8,9) are used. Acceptable type specific classification is not possible using whole virus lysate, the commonly used antigen of early HSV serology kits. The most commonly used type specific protein is glycoprotein G. **ImmunoFLOW™ HSV uses HSV gG type 2 recombinant protein.**

### HSV IgG ImmunoFLOW™

T ImmunoFLOW is an immunoassay consisting of a cassette and three reagents. The cassette contains a paper matrix (for example, nitrocellulose) and an absorbent material. The paper matrix was manufactured with three "dots", each contain an antigen (e.g., positive control, analyte 1, analyte 2). A body fluid (e.g., serum) is applied to the triangular opening and allowed to flow through the paper matrix into the absorbent. To assure assay specificity, a wash reagent is applied and flows into the absorbent material. Finally, gold particles attached to an immunological reagent is applied and absorbed. If specific antibody binding occurs, immunologically active gold particles will bind and cause a red/pink color formation.

### PERFORMANCE:

#### HSV IgG ImmunoFLOW™

Rapid Test color immunochromatographic, capillary flow technology

Sensitivity: 100%

Relative Specificity: 100%

Reproducibility: 100%