fr HYDATIDOSE FUMOUCZ®

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en ECHINOCOCCOSIS FUMOUCZ®

INTENDED USE: SERODIAGNOSIS OF ECHINOCOCCOSIS BY INDIRECT HAEMAGGLUTINATION

ECHINOCOCCOSIS FUMOUCZ® is an indirect haemagglutination test for the quantitative detection of anti-*Echinococcus granulosus* antibodies in serum samples.

PRINCIPE :
ECHINOCOCCOSIS FUMOUCZ® principle is based on indirect haemagglutination. Sensitized red blood cells are composed of sheep red blood cells coated with an *Echinococcus granulosus* antigen. Serum antibodies against *Echinococcus granulosus* are revealed by an agglutination of the sensitized red blood cells: a reddish-brown film can be observed in the well. In the absence of specific antibodies, these red blood cells deposit, forming a ring in well bottom. The unsensitized red blood cells ensure the reaction specificity and allow the elimination of interferences due to natural anti-sheep agglutinins (Forssman heteroantibodies, infectious mononucleosis antibodies...). The reaction is performed in U-microplate. The test procedure is easy and rapid. The results are obtained in 2 hours.

KIT CONTENT:

- R1**: 1 vial of 2.2 mL of sensitized red blood cells
- BUF**: 1 vial of 55 mL of phosphate buffer pH 7.2
- CONTROL +**: 1 vial of 0.2 mL of titrated positive control
- MICROPLATE**: 2 microplates with a U-bottom
- R2**: 1 vial of 1 mL of non-sensitized red blood cells
- R3**: 1 vial of 2 mL of adsorbent
- CONTROL -**: 1 vial of 0.2 mL of negative control
- DROPPER**: 2 special droppers

Notes : - Use exclusively the droppers provided in the kit.
- Do not interchange droppers between **R1** and **R2** reagents.

MATERIAL REQUIRED BUT NOT PROVIDED:

- Multi-channel micropipette delivering 50 µL
- Micropipette delivering 50 µL, micropipette 100-1000 µL
- Container for contaminated wastes
- Haemolysis tubes
- Centrifuge

STORAGE CONDITIONS:

Reagents are ready-to-use. Store at +2°...+8°C, until the expiry date indicated on the box. Do not freeze.

SAMPLES COLLECTION / HANDLING / STORAGE:

Use fresh withdrawn serum. The serum samples should be stored at +2°...+8°C for a period of 72 hours. If the sera are not tested within 72 hours after their collection, they must be stored (-20°C). It is recommended to prepare aliquots to avoid repeated freezings and defrostings. Do not perform a heat pre-treatment of serum. Do not use haemolysed, turbid or contaminated serum.

WARNINGS AND PRECAUTIONS:

- For in vitro diagnostic use and only for professional use. Tests are for single use only.
- Reagents and controls from different lots must not be interchanged.
- Follow the instructions for use. Avoid contact of reagent with skin, eyes and mucous membranes. Do not ingest.
- Allow the reagents and samples to return to room temperature prior to testing.
- Carefully shake the red blood cells suspensions before use.
- When delivering **R1** and **R2** reagents, make sure that the dropper needle is perfectly vertical. Check the absence of air bubbles in the drops to ensure constant delivery volumes.
- In case of accidental spill of reagents, clean the surface with absorbent paper and water. In case of spill of serum, clean the surface with absorbent paper and bleach.
- The samples, reagents as well as the contaminated materials and products must be eliminated in a container for contaminated wastes, according to the prevailing recommendations and regulations.
- Controls contain < 0.1 % sodium azide which may react with lead and copper plumbing to form explosive metal azides. Azide built-up can be avoided by flushing with large volumes of water following the disposal reagents.
- All reagents, except **BUF** reagent, contain material of animal origin. As a consequence, they have to be handled with care as hazardous components.

TEST PROCEDURES:

Before use, allow reagents and samples to return to room temperature.

- Preparation of 1:40 stock dilution of test serum**
Deliver in a haemolysis tube and mix: 50 µL of test serum and 1.95 mL of **BUF** reagent.
 - Test execution on microplate**
 - By means of a multi-channel micropipette, deliver 50 µL of **BUF** reagent in 8 wells of the microplate.
 - Add 50 µL of serum stock dilution in the 1st well. Mix with **BUF** reagent and transfer by means of a micropipette 50 µL from the 1st well into the 2nd well, from the 2nd into the 3rd, and so on until the 6th well. Discard 50 µL from the 6th well. The following dilutions are obtained:
- | Well number (see annex) | Well I | Well II | Well III | Well IV | Well V | Well VI |
|-------------------------|--------|---------|----------|---------|--------|---------|
| Dilution | 1:80 | 1:160 | 1:320 | 1:640 | 1:1280 | 1:2560 |
- Add 50 µL of serum stock dilution into the 7th well. Mix with **BUF** reagent and discard 50 µL. This 1:80 dilution constitutes the serum control. It serves for the detection of natural anti-sheep agglutinins, which may occur in certain sera.
 - Carefully shake **R1** and **R2** reagents.
 - Distribute 1 drop of **R1** reagent in the first 6 wells.
 - Distribute 1 drop of **R2** reagent in the 7th well (serum control - well VII - see annex).
 - Distribute 1 drop of **R1** reagent in the 8th well (reagent control - well VIII - see annex). Its role is to control the validity of **BUF** reagent and **R1** reagent.
- Note : Set up only one reagent control per assays series.
- Very carefully homogenize wells content manually, by lateral thrumming on the edges of the microplate, placed flatwise or with a vibrator-shaker for microtitration plates. Do not use orbital shaker. Then allow the plate to remain motionless, protected from vibrations.
 - Read the reaction 2 hours later.

ADSORPTION OF NATURAL ANTI-SHEEP AGGLUTININS IN CASE OF AGGLUTINATION OF SERUM CONTROL:

- Deter in a tube and mix: 0.1 mL of serum and 0.3 mL of **R3** reagent .
- Incubate for 60 min at room temperature.
- Centrifuge at 800 g for 15 min.
- Collect the supernatant fluid; the serum is then 1:4 diluted.
- Dilute the supernatant fluid 1:10 with the **BUF** reagent to obtain a 1:40 adsorbed stock dilution.
- Follow the protocol of "Test execution on microplate" by replacing the stock dilution by the adsorbed stock dilution.

INTERNAL QUALITY CONTROL:

- Each **ECHINOCOCCOSIS FUMOUCZ®** kit includes **CONTROL -** and **ltered CONTROL+**. They are ready-to-use and have to be run as the samples. They validate the test. The **CONTROL+** titer must be equal to ± one dilution from the one mentioned on the vial label. No haemagglutination must be observed with the **CONTROL -**. If not, the test is not valid.
- The serum control must give a negative reaction (ring). In case of haemagglutination of this control, it is necessary to repeat the test after elimination of natural anti-sheep agglutinins by adsorption.
- The reagent control must give a negative reaction (ring). In case of haemagglutination of this control, **ECHINOCOCCOSIS FUMOUCZ®** reagent could not be used.

READING OF RESULTS:	NO HAEMAGGLUTINATION Presence of a more or less wide ring in well bottom.	NEGATIVE REACTION
	HAEMAGGLUTINATION Presence of a reddish-brown film in well bottom; sometimes, presence of a thin peripheral ring.	POSITIVE REACTION

Haemagglutination picture, obtained with a 1:1280 positive serum, is presented in annex.

The titer is given by the first dilution showing a wide peripheral ring.

- Notes :
- In case of positive reaction in the first 6 wells, make higher dilutions to determine the limit haemagglutination titer.
 - Some sera, with very high antibodies titer, may give zone phenomenon (with film retraction) at the first dilutions, which disappears at increasing dilutions.
 - The quality of the reagents allows the performance of the reaction in the evening, with reading the next morning, providing the microplate remains motionless and protected from vibrations.

INTERPRETATION OF THE RESULTS:

TITER < 1:160	Non significant reaction. Probable absence of echinococcosis. Repeat the test 2 or 3 weeks later and associate a counter-immunoelectrophoresis or an immunoelectrophoresis.
TITER = 1:160	Equivocal reaction. Repeat the test 2 or 3 weeks later and associate a counter-immunoelectrophoresis or an immunoelectrophoresis.
TITER ≥ 1:320	Significant reaction in favour of acute echinococcosis.

PERFORMANCE CHARACTERISTICS:

Echinococcosis disease development is often slow and not very symptomatic. Therefore, the sero-diagnosis is of high interest. The **ECHINOCOCCOSIS FUMOUCZ®** reagent is composed of red blood cells sensitized with an *Echinococcus granulosus* antigen. It ensures sensibility and specificity to the reaction. A study of 221 human sera has shown a sensitivity of 93.0 % (whichever cyst location) and a specificity of 94.9 %. Comparative study with IFI and with ELISA test has shown the complementarity of these various test (2). In any case, diagnosis should be made using the results of this test together with the other clinical, epidemiological and laboratory findings.

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