SARS-CoV-2 Triplex PCR kit

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1. INTENDED USE

SARS-CoV-2 Triplex PCR kit is one-step reverse transcription real-time PCR multiplex test for detection of SARS-CoV-2 nucleic acids from nasopharyngeal swabs, oropharyngeal swabs and saliva samples from individuals suspected of COVID-19.

2. PRINCIPLE OF THE TEST

The test is based on the principle of one-tube subsequent reactions of reverse transcription and real-time PCR with fluorescent probe detection.

During first reaction, purified viral RNA is reverse transcribed into cDNA. During second reaction, target sequences are amplified from cDNA template with target-specific primers. During amplification reaction, target-specific fluorescent probes are bound to amplified fragments and then degrade by polymerase 5' → 3' exonuclease activity. Released fluorophores generate signal that is detected and analysed by real-time PCR thermocycler.

SARS-CoV-2 Triplex PCR kit simultaneously detects three targets, in ORF1ab, ORF8 and N protein coding regions, each on separate detection channel. This allows increasing the robustness of the test in case of possible arising mutations in the regions where primers or probes bind.

3. KIT FORMATS AND CONTENTS

3.1. Form S (PCR-Mix SC2 reaction mix under paraffin in 0.1 mL low profile 8-tube strips)

Label	Component	REF 89-01		REF 89-01-24	
		Volume; μL	Quantity	Volume; μL	Quantity
SC2-Mix-S	PCR-Mix SC2, Form S	10	13×8×0.1 mL tube	10	3×8×0.1 mL tubes
Taq+RT	Taq DNA-polymerase with Reverse Transcriptase	510	1 tube	125	1 tube
PC	Positive Control	Lyophilized powder	1 tube	Lyophilized powder	1 tube
NC	Negative Control	150	1 tube	150	1 tube

The Form S of the product was validated to be used with the following real-time PCR cyclers: Bio-Rad CFX96 (also can be used with CFX96 Touch), DNA-Technology DTprime (corresponds to DTlite, SaCycler), and Roche Light Cycler 96. While using other real-time PCR cyclers compatible with 0.1 mL low profile PCR tubes, contact the kit manufacturer for more information concerning compatibility.

3.2. Form T (PCR-Mix SC2 reaction mix under paraffin in individual 0.2 mL regular profile tubes)

Label	Component	REF 89-02		REF 89-02-24	
		Volume; μL	Quantity	Volume; μL	Quantity
SC2-Mix-T	PCR-Mix SC2, Form T	10	100 x 0.2 mL tubes	10	24 x 0.2 mL tubes
Taq+RT	Taq DNA-polymerase with Reverse	510	1 tube	125	1 tube

	Transcriptase				
PC	Positive Control	Lyophilized powder	1 tube	Lyophilized powder	1 tube
NC	Negative Control	150	1 tube	150	1 tube

The Form T of the product was validated with following real-time PCR cyclers: Qiagen Rotor-Gene Q (also can be used with Corbett Research Rotor-Gene 3000/6000), Bio-Rad CFX96 (also can be used with CFX96 Touch), DNA-Technology DTprime (corresponds to DTlite, SaCycler). While using other real-time PCR cyclers compatible with regular profile PCR tubes, contact the kit manufacturer for more information concerning compatibility.

3.3. Form F (liquid PCR-Mix SC2 reaction mix is in 2 mL tube)

Label	Component	REF 89-03		REF 89-03-24	
		Volume; μ L	Quantity	Volume; μ L	Quantity
SC2-Mix-F	PCR-Mix SC2, Form F	1020	1 tube	250	1 tube
Taq+RT	Taq DNA-polymerase with Reverse Transcriptase	510	1 tube	125	1 tube
PC	Positive Control	Lyophilized powder	1 tube	Lyophilized powder	1 tube
NC	Negative Control	150	1 tube	150	1 tube

The Form F of the product was validated to be used with following real-time PCR cyclers: Qiagen Rotor-Gene Q (also can be used with Corbett Research Rotor-Gene 3000/6000), Bio-Rad CFX96 (also can be used with CFX96 Touch), DNA-Technology DTprime (corresponds to DTlite, SaCycler), Roche Light Cycler

96. While using other real-time PCR cyclers contact the manufacturer for more information concerning compatibility.

4. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Biosafety cabinet;
- Real-Time PCR cycler, dependent of the kit form used, e.g. Bio-Rad CFX-96 , DNA-Technology DTprime (corresponds to DTlite, SaCycler), Qiagen Rotor-Gene Q (Corbett Research Rotor-Gene 6000), Roche Light Cycler96;
- Vortex mixer;
- Microcentrifuge for 1.5-2 mL eppendorf tubes;
- PCR tubes/microwell plates, if using kit form F;
- Microcentrifuge for 0.1-0.2 mL PCR tubes, 8-tube strips or microwell plates, depending on kit form used;
- Two separate pipette sets for DNA/RNA-free and DNA/RNA-containing components;
- Filter tips (10, 100; 200 and 1000 µL);
- Disposable safe-lock 1.5-2 mL microcentrifuge tubes, DNase/RNase free;
- Magnetic DNA/RNA extraction kit (Ref 80-06) or Prime DNA/RNA Extraction kit (Ref 80-07);
- Separate working areas for DNA/RNA extraction (Sample Preparation Area), PCR preparation (Reagent Preparation Area) and for real-time PCR amplification (Amplification Area). Each area must have own set of pipettes, instruments and laboratory clothing.

5. STORAGE CONDITIONS AND STABILITY OF THE KIT

SARS-CoV-2 Triplex PCR kit should be stored at -18...-30 °C, preferably in the original kit box, until the expiration date. The manufacturer guarantees the functionality of the kit components up to 15 freezing and thawing cycles.

The expiry date of the kit is stated on the kit box label, expiry date for each component is indicated on the respective tube label. Storage at +25 °C is allowed for no more than 5 days (any kit format).

6. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE

Total DNA/RNA can be extracted from nasopharyngeal swabs; oropharyngeal swabs and saliva, collected for transportation and storage in recommended transport medium. Please, ensure that sample collection, transportation and storage conditions meet all necessary local regulations.

Type of specimen	Recommendations for collection and transportation	Recommendations for storage
Nasopharyngeal and oropharyngeal swabs, saliva	Specimens are collected with recommended disposable collection devices and transported in 0.5 mL of 0.9 % NaCl solution (saline) or other suitable transport medium refrigerated at +2...+8 °C	<ul style="list-style-type: none"> ○ At 18-25 °C for ≤ 48 hours ○ Refrigerated at 2-8 °C for ≤ 4 days ○ Frozen at -18...-35 °C for > 4 days

7. ASSAY PROCEDURE

The complete procedure consists of three stages:

I. RNA extraction

- II. Real-Time RT-PCR (one-step)
- III. Data analysis and interpretation

I. RNA extraction

NOTE: All manipulations with samples and DNA/RNA extraction must be done in Sample Preparation Area.

NOTE: The quality of the extracted RNA has a profound impact on the performance of the entire test system. It has to be ensured that the system used for RNA extraction yields RNA of quality and quantity suitable for real-time PCR. SARS-CoV-2 Triplex PCR kit was validated with the following kits and systems for DNA/RNA extraction, and are strongly recommended by the manufacturer:

- Magnetic DNA/RNA extraction kit (Astra Biotech, Ref 80-06)
- Prime DNA/RNA Extraction kit (Astra Biotech, Ref 80-07)

The compatibility of other DNA/RNA extraction procedures for use with SARS-CoV-2 Triplex PCR kit has to be validated by the user.

- 1) Prepare and label required quantity of 1.5-2 mL microcentrifuge tubes equal to the number of samples and 1 tube for negative extraction control labeled as “NEC”.

NOTE: Each extraction procedure should be performed using “NEC”.

NOTE: You do not need to put in each sample Internal Control **IC**, because test system uses endogeneous internal control (transcript of one of human «housekeeping genes»).

2) Add 100 µL **NC** into the microcentrifuge tube marked “NEC”.

NOTE: The Elution buffer or Negative control contained in DNA/RNA extraction kit can be used instead of negative extraction control.

3) Follow the manual of DNA/RNA extraction kit used.

II. Real-Time RT-PCR

IMPORTANT: Positive and Negative Controls must be included in each run.

Form S/T

Kit components preparation

NOTE: Kit components preparation has to be done in the Reagent Preparation Area

- 1) Take out from the freezer and thaw the following components of the kit:
 - **PC** (positive control);
 - **NC** (negative control);
 - **Taq+RT** (Taq polymerase with Reverse Transcriptase);
 - Required quantity of 8-tube strips or individual PCR tubes with **SC2-Mix-S/T** for tested samples, plus 1 tube for “NEC” and 1 tube for “PC”.

The rest of components put back in the freezer.

2) Transfer **PC** and **NC** to the Sample Preparation Area.

PCR mix preparation

NOTE: PCR mix preparation has to be done in the Reagent Preparation Area. Do not prepare the PCR mix in the same area as the DNA/RNA extraction or PCR amplification to prevent contamination.

- 1) Label prepared PCR tubes with **SC2-Mix-S/T**.

NOTE: When using real-time PCR cyclers with fluorescent detection through the tube cap, label tubes on the edge of the cap or on side wall. When using real-time PCR cyclers with fluorescent detection through the side wall, label tubes on the cap.

- 2) Gently vortex **Taq+RT**, avoiding foaming, and centrifuge at 1500-2400 rpm for 2-3 s.
- 3) Place 10 µl of **Taq+RT** into each tube avoiding tip contact with paraffin.

Dispensing of test samples and controls

NOTE: Dispensing of test samples and controls has to be performed in the Sample Preparation Area.

NOTE: Use a separate filter tip for each DNA sample.

- 1) Centrifuge the RNA extracted from the specimens and “NEC” (**NC** extract) at 13000 rpm for 1 min before adding them to the PCR tubes.
- 2) Add 10 µl of specimen RNA and “NEC” to corresponding PCR tubes avoiding tip contact with paraffin.
- 3) Add 50 µL of **NC** to the **PC** tube.
- 4) Vortex and centrifuge the **PC** at 1500-2400 rpm for 2-3 s.
- 5) Add 10 µL of **PC** to the PCR tube labeled “PC” avoiding tip contact with paraffin.

6) Close PCR tubes and centrifuge at 1500-2400 rpm for 2-3 s.

Real-Time PCR Amplification

NOTE: *Real-Time PCR Amplification has to be performed in the Amplification area. Separate room for this area is recommended. Follow the manufacturer instruction before using PCR Real-Time cycler.*

- Load PCR tubes into a real-time PCR cycler.
- Using cycler's software, fill in the location of the samples. Set up fluorescent detection on the channels ROX/Orange, FAM/Green, HEX/Yellow and Cy5/Red.
- Run the PCR protocol.

Stage	Temperature, °C	Time	Data collection	Number of cycles
Hold	50	15 min	-	1
Hold	94	3 min	-	1
Cycling	94	10 s	-	5
	60	20 s	-	
Cycling	94	10 s	-	45
	60	20 s	ROX/ Orange FAM/ Green HEX/ Yellow Cy5/ Red	

Form F

Kit components preparation

NOTE: *Kit components preparation has to be done in the Reagent Preparation Area*

- 1) Take out from the freezer and thaw the following components of the kit:

- **PC** (positive control);
 - **NC** (negative control);
 - **Taq+RT** (Taq polymerase with Reverse Transcriptase);
 - **SC2-Mix-F** (PCR buffer with dNTPs, primers and probes).
- 2) Transfer **PC** and **NC** to the Sample Preparation Area.

PCR mix preparation

NOTE: *PCR mix preparation has to be performed in the Reagent Preparation Area. Do not prepare the PCR mix in the same area as the DNA/RNA extraction or PCR amplification to prevent contamination.*

- 1) Prepare and label required quantity of 0.1-0.2 mL PCR tubes/strips/microwell plates for tested samples, plus 1 tube/well for “NEC” and 1 tube/well for “PC”.

Note: *when using real-time PCR cyclers with fluorescent detection through the tube cap, label tubes on the edge of the cap or on side wall. When using real-time PCR cyclers with fluorescent detection through the side wall, label tubes on the cap.*

- 2) Vortex and centrifuge the **SC2-Mix-F** at 1500-2400 rpm for 2-3 s.
- 3) Gently vortex **Taq+RT**, avoiding foaming, and centrifuge at 1500-2400 rpm for 2-3 s.
- 4) To prepare Master Mix, in 1.5-2 mL eppendorf tube, mix $10 \times N$ μ L of **SC2-Mix-F** and $5 \times N$ μ L of **Taq+RT**, where N is the total number of reactions.

For example, to run 10 samples, and «PC» and «NEC» (12 reactions total) you should mix 120 μ L of **SC2-Mix-F** and 60 μ L of **Taq+RT**.

NOTE: Tube of **SC2-Mix-F** is sufficient for 100 reactions plus 2 additional reactions to account for pipetting errors. To set up exactly 100 reactions, add 510 μ L of **Taq+RT** directly into unused tube of **SC2-Mix-F**.

5) Vortex prepared Master Mix and spin briefly to remove drops from the inside of the lid.

NOTE: Do not store Master Mix (**SC2-Mix-F** plus **Taq+RT**) for more than 1 hour at +2...+8 °C.

6) Put 15 μ L of Master Mix into each PCR tube/ well.

Dispensing of test samples and controls

NOTE: Dispensing of test samples and controls have to be performed in the Sample Preparation Area.

NOTE: Use a separate filter tip for each RNA sample.

7) Centrifuge the RNA extracted from the specimens and “NEC” (**NC** extract) at 13000 rpm for 1 min before adding them to the PCR tubes.

8) Add 10 μ L of specimens RNA and “NEC” in corresponding PCR tubes/wells.

9) Add 50 μ L of **NC** to the **PC** tube.

10) Vortex and centrifuge the **PC** at 1500-2400 rpm for 2-3 s.

11) Add 10 μ L of **PC** to the PCR tube/well labeled “PC”.

12) Close the tubes and centrifuge at 1500-2400 rpm for 2-3 s.

Real-Time PCR Amplification:

NOTE: Real-Time PCR Amplification has to be performed in the Amplification area. Separate room for this area is recommended. Follow the manufacturer’s instructions before using the PCR Real-Time cycler.

- Load PCR tubes into a real-time PCR cycler.
- Using cycler's software, fill in the location of the samples. Set up fluorescent detection on the channels ROX/Orange, FAM/Green, HEX/Yellow and Cy5/Red.
- Run the PCR protocol.

Stage	Temperature, °C	Time	Data collection	Number of cycles
Hold	50	15 min	-	1
Hold	94	3 min	-	1
Cycling	94	10 s	-	5
	60	20 s	-	
Cycling	94	10 s	-	45
	60	20 s	ROX/ Orange FAM/ Green HEX/ Yellow Cy5/ Red	

III. Data Analysis

1) General conditions:

Fluorescent dyes ROX, HEX and Cy5 are used to detect the SC2 RNA and FAM is used to detect the **IC** RNA.

Analysis of the results can be performed with Real-Time PCR thermal cycler's software according to its instruction for use. For the validated cyclers detailed instructions concerning protocol set up and results analysis are given in «Astra Biotech: Real-Time PCR kits Guidelines».

Analysis of the results can be performed only by using following parameters: Cq **ORF1**, Cq **ORF8**, Cq **N**, Cq **PC** and Cq **IC** that are specified for each lot in the Quality Control Sheet provided with every kit.

Cq **ORF1** = critical value of quantification cycle for SC2 in the ROX channel, is crucial for detection positive results;

Cq **ORF8** = critical value of quantification cycle for SC2 in the HEX channel, is crucial for detection positive results;

Cq **N** = critical value of quantification cycle for SC2 in the Cy5 channel, is crucial for detection positive results;

Cq **PC** = critical value of quantification cycle for Positive Control (PC) in the ROX, HEX and Cy5 channels is used to determine whether kit components for reverse transcription and PCR amplification and detection work properly is crucial for detection negative results;

Cq **IC** = critical value of quantification cycle for Internal Control (IC) in the FAM channel is used to control all stages of the analysis, including quality of RNA extraction procedure.

The results of test run are valid if:

- Cq PC \leq Cq **PC** in the ROX, HEX and Cy5 channels. In the case when Cq value for PC is $>$ Cq **PC**, then all the specimens defined as positive in this run should be considered as valid, whereas all the specimens defined as negative should be considered as invalid.
- Cq of NEC $>$ Cq **ORF1** in the ROX channel and Cq of NEC $>$ Cq **IC** in the FAM channel, but there are some exceptions (for more information see 2) Data interpretation).

- Cq of NEC $>$ Cq **ORF8** in the HEX channel and Cq of NEC $>$ Cq **IC** in the FAM channel, but there are some exceptions (for more information see 2) Data interpretation).
- Cq of NEC $>$ Cq **N** in the Cy5 channel and Cq of NEC $>$ Cq **IC** in the FAM channel, but there are some exceptions (for more information see 2) Data interpretation).

2) Data interpretation

Specimen is considered as positive for RNA of SC2 if for any two of three targets:

$Cq \leq Cq$ **ORF1** in the ROX channel

$Cq \leq Cq$ **ORF8** in the HEX channel

$Cq \leq Cq$ **N** in the Cy5 channel

Cq of PC $\leq Cq$ **PC** in the ROX channel

even if $Cq > Cq$ **IC** in the FAM channel or is not determined.

Specimen is considered as negative for RNA of SC2:

$Cq > Cq$ **ORF1** in the ROX channel or is not determined

$Cq > Cq$ **ORF8** in the HEX channel or is not determined

$Cq > Cq$ **N** in the Cy5 channel or is not determined

$Cq \leq Cq$ **IC** in the FAM channel,

$Cq \leq Cq$ **PC** in the ROX channel.

Specimen is considered as invalid if for any two of three targets:

$Cq > Cq$ **ORF1** in the ROX channel or is not determined

$Cq > Cq$ **ORF8** in the HEX channel or is not determined

$Cq > Cq$ **N** in the Cy5 channel or is not determined

but any of the following is true:

$Cq > Cq_{IC}$ in the FAM channel or is not determined

$Cq > Cq_{PC}$ in the ROX, HEX, Cy5 channel or is not determined

In case of contamination:

$Cq_{NEC} \leq Cq_{ORF1}$ in the ROX channel and/or

$Cq_{NEC} \leq Cq_{ORF8}$ in the HEX channel and/or

$Cq_{NEC} \leq Cq_N$ in the Cy5 channel. Measures must be taken to eradicate the source of contamination. In case of contamination, all positive specimens satisfying condition $Cq \geq (Cq_{NEC} - 5)$ in the ROX, HEX, Cy5 channels are considered as invalid, the other results are considered as valid.

For all specimens with invalid results it is recommended to perform the analysis from the RNA extraction stage.

8. PERFORMANCE CHARACTERISTICS OF THE KIT

Analytical sensitivity:

The detection limit of the kit evaluated via probit analysis is 1000 copies of SC2 genomic RNA per 1 mL of sample (or ca. 10 copies per reaction). This limit is defined for confidence level 95 %. It is assumed that 100 µl of specimen was taken for DNA/RNA extraction and elution volume was 60 µl.

Analytical specificity:

In silico analysis of primer and probe sequences showed no homology with genome sequences of the following organisms: *Humancoronavirus 229E*, *Human coronavirus OC43*, *Human*

coronavirus HKU1, Human coronavirus NL63, SARS-coronavirus, MERS-coronavirus, Adenovirus (e.g. C1 Ad. 71), Human Metapneumovirus (hMPV), Parainfluenzavirus 1-4, Influenza A & B, Enterovirus (e.g. EV68), Respiratory syncytial virus, Rhinovirus, Chlamydia pneumoniae, Haemophilus influenzae, Legionella pneumophila, Mycobacterium tuberculosis, Streptococcus pneumoniae, Streptococcus pyogenes, Bordetellapertussis, Mycoplasma pneumoniae, Pneumocystis jirovecii (PJP), Candida albicans, Pseudomonas aeruginosa, Staphylococcus epidermis, Staphylococcus salivariuss.

Diagnostic sensitivity and specificity:

Diagnostic evaluation of SARS-CoV-2 Triplex PCR kit was performed on 210 clinical samples (nasopharyngeal swabs, oropharyngeal swabs, and saliva). RNA from the samples was extracted with Prime DNA/RNA Extraction kit.

Results of SARS-CoV-2 Triplex PCR kit were compared with those of reference diagnostic kit for SARS-CoV-2 detection (Reference kit) and are presented in the table below:

SARS-CoV-2 Triplex PCR kit results	Reference kit Positive (47)	Reference kit Negative (163)
Positive	47	0
Negative	0	163

Diagnostic sensitivity is 100%.

Diagnostic specificity is 100%.

9. LIMITATIONS OF THE METHOD

Any clinical diagnosis should not be based on the results of *in vitro* diagnostic methods alone. For final diagnosis a physician is supposed to consider all available clinical and laboratory tests findings.

10. SAFETY PRECAUTIONS

Taq DNA-polymerase+RT **Taq+RT**, Positive Control **PC** are harmful, irritant and can cause serious eye damage; the following precautions should be observed:

- Wash the hands thoroughly after handling.
- Do not eat, drink or smoke when using this product.
- Wear protective gloves/ protective clothing/eye protection/ face protection.
- Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.

- Use only outdoors or in a well-ventilated area.
 - Store in a well-ventilated area.
 - IF INHALED: remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a poison center or doctor/physician if you feel unwell. Do not induce vomiting.
 - IF SWALLOWED: call a poison center or doctor/physician if you feel unwell. Rinse mouth.
 - IF ON SKIN: remove/ take off immediately all contaminated clothes, rinse skin with plenty of soap and water/ shower. Call a poison center or doctor/physician if you feel unwell.
 - IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.
 - Wash contaminate clothing before reuse.
-
- **This kit is for in vitro diagnostic use only.** The operator should thoroughly follow the manual to obtain reliable data. This instruction manual is valid only for the present kit with the listed contents. Any exchange of the kit components is not allowed by CE regulations.
 - Test should only be performed by skilled personnel who follow GLP (Good Laboratory Practice) guidelines.
 - Don't use the kit after its expiration date.
 - Do not pool reagents from different lots or from different tubes of the same lot. After use, immediately close all tubes in order to avoid leakage.
 - We recommend setting up blank protocols of PCR before the analysis. We recommend filling all necessary sample information in PCR cycler software beforehand.

- PCR technology is extremely sensitive and prone to contamination issues. Therefore, set up three separate working areas for a) sample preparation b) PCR reagent preparation and c) PCR amplification. For each working area an own set of pipettes, instruments and protective clothing should be reserved.
- Use sterile filter tips for pipetting and use special PCR pipettes for aerosol free pipetting.
- Always use fresh tips to pick different solutions.
- Routinely decontaminate your pipettes and the laboratory benches.
- After usage, all the reagents and test components should be discarded according to local regulations.

Taq DNA-polymerase+RT **Taq+RT**, Positive Control **PC** are harmful for aquatic life. Avoid release to the environment.

-  Any materials of animal origin used for manufacturing of kit components were found safe for humans. However, none of known laboratory test guarantees total absence of the hazardous agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage, all kit components, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving or disinfectant treatment. After sterilization all kit

components and disposable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.

-  The following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves while performing the assay;
 - never pipette material by mouth;
 - be careful while handling **PC**, avoid its splashing;
 - in case of any reagent or sample spillage, wipe up the spills promptly and wash affected area thoroughly using suitable decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

DISTRIBUITO IN ITALIA DA

MERIDIAN HEALTHCARE srl
Via Caronda, 446 SC/A - 95100 CATANIA - ITALIA
Tel.: +39 095 7256869 Fax: +39 095 7254454
e-mail: info@meridianhealthcare.it
web: www.meridianhealthcare.it

 Meridian Healthcare®