BIOSYNEX

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AMNIOQUICK[®] CARD IGFBP-1

REF. 1090005

BIOSYNEX

12, rue Ettore Bugatti – CS 28006

67038 STRASBOURG Cedex – France

Tél : +33 3 88 77 57 00 - Fax : +33 3 59 81 21 74

Site : www.biosynex.com - Mail : info@biosynex.com

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MERIDIAN HEALTHCARE srl

Via G. Guglielmino, 68

95030 Tremestieri Etneo (Catania) Italy

Tel. +39 095 7256869 Fax: +39 095 7254454

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


BIBLIOGRAFIA

1. Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550.
2. Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851.
3. Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610.
4. Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418.
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915.
6. Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300.
7. Bell SC, Keyte JW; N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein – Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988;1202-1204.
8. Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81.
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic 2006;73: 659-664.
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150.
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'Insulin-like growth factor-binding protein-1 et de la fibronectine foetale dans les sécrétions cervicovaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999 ;36 :388-390.
12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-I; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161.
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 :12 :233-240.
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005;121:159-163.

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