Auro•Dex®

Helicobacter Pylori MultiTest

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BIBLIOGRAFIA

1. NCCLS Doc M29-P 1988. Protection of laboratory workers from infectious disease transmitted by blood and tissue, Proposed guideline. National Committee for Clinical Laboratory Standards. Villanova, PA.
2. CDC/NIH Manual. 1988. Biosafety in Microbiological and Biomedical Laboratories, 2nd ed. pp 12-16
3. Jones, D.M., A.M. Lessells, J. Eldridge: Campylobacter-like organisms on the gastric mucosa culture, histological and serological studies. J. Clin. Pathol., 37, 1002, 1984.
4. Buck, G.E., W.K. Gourley, K. Subramanyam, J.M. Latimer, A.R. DiNuzzo: Relation of Campylobacter pyloridis to gastritis and peptic ulcer. J. Infect. Dis., 153, 664. 1986.
5. Marshall, B.J., H. Royce, D.I. Annear, C.S. Goodwin, J.W. Pearman, J.R. Warren and J.A. Armstrong. 1984. Original Isolation of Campylobacter pyloridis from human gastric mucosa. Microbics Lett. 25:83-88.
6. Marshall, B.J., and J.R. Warren. 1984. Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. Lancet. i:1311-1314.
7. Blaser, M.J. Helicobacter pylori and the pathogenesis of gastroduodenal inflammation. Journal of Infectious Diseases. 1990. 161:626-633.
8. Graham, D.Y. 1989. Campylobacter pylori and peptic ulcer disease. Gastroenterology. 96:615-625.
9. Dooley C.P., P.L. Fitzgibbons, H. Cohen, M.D. Appleman, G.I. Perez-Perez and M.J. Blaser. 1989. Prevalence of Helicobacter pylori infection and histologic gastritis in asymptomatic persons. N. Eng. J. Med. 321:1562-66.
10. Booth L., Holdstock G, MacBride H, et al. Clinical importance of Campylobacter pyloridis and associated serum IgG and IgA antibody responses in patients undergoing upper gastrointestinal endoscopy. J Clin Pathol 1986;39:215-9.
11. Yamamoto, I., Y. Fukuda, T. Mizuta, M. Fukada, T. Nishigami, and T. Shimoyama. 1995. Serum anti-Helicobacter antibodies and gastritis. J. Clin. Gastroenterol. 21:5164-8.
12. Andersen, L., S.J. Rosenstock, O. Bonnevie, and T. Jorgensen. 1996. Seroprevalence of immunoglobulin G, M, and A antibodies to Helicobacter pylori in an unselected Danish population. Am J. Epidemiol. 143(11):1157-64.
13. Tsavaris N., et al. 1998. Preliminary evaluation of the potential prognostic value of serum levels of immunoglobulins (IgA, IgM, IgG, IgE) in patients with gastric ulcer. Int. J. Biol. Markers. 13(2):87-91.
14. Aromaa, A., et al. 1996. Circulating anti-Helicobacter pylori immunoglobulin A antibodies and low serum pepsinogen I level are associated with increased risk of cancer. Am. J. Epidemiol. 144(2):142-9.
15. Kosunen, T.U., K. Seppala, S. Sarna, and P. Sipponen. 1992. Diagnostic value of decreasing IgG, IgA, and IgM antibody titres after eradication of Helicobacter pylori. Lancet. 339(8798):893-5.
16. Andersen, L.P., A.A. Wewer, K.M. Christiansen, M. Tvede, J.P. Hansen, F.W. Henriksen, and P.A. Krasilnikoff. 1994. The humoral immune response to Helicobacter pylori infection in children with recurrent abdominal pain. APMS. 102(6):457-64.
17. Borody T.J., Shortis NP, Reyes, E. Eradication therapies for Helicobacter pylori. J. Gastroenterol 1998;33 Suppl 10:53-6.
18. Okada M., Oki K, Shirotani T, Seo M, Okabe N, Maeda K, Nishimura H, Ohkuma K, Oda K. A new quadruple therapy for the eradication of Helicobacter pylori. Effect of pretreatment with omeprazole on the cure rate. J. Gastroenterol 1998 Oct;33 (5):640-5.

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EUROPE CONTACT
Meridian Healthcare[®]

Meridian Healthcare srl
68, via G. Guglielmino
Tremestieri Etneo, CT 95030
Phone: +39 095 725 68 69
FAX: +39 095 725 44 54

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Awareness Technology Inc.