MICROALBUMINA HSA TEST RAPIDO
REF.5001C31-5 25 test

USO PREVISTO

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BIBLIOGRAFIA

1. Borch, Johnsen K, Wenzel, H Viberti, G. C, and Mogensen C. E, Is screening and intervention for microalbuminuria worthwhile in patients with insulin dependent diabetes Br, Med. J 306, 1722-1723, 1993.
2. Whiteside, C. I Clin. Biochem. Detection of progressive diabetic nephropathy role of microalbuminuria determination. 26, 311-313, 1993.
3. Mogensen, C. E. Proceedings of a meeting on epidemiology of microalbuminuria, J Diabet Complications 8, 133-192, 1992.
4. De Cotret PR. Relationships among diabetes, microalbuminuria, and ACE inhibition. J Cardiovasc

Pharmaco1, 32(Supp1.2), S9-S17), 1998.

5. Mogensen Ce, et al, Prevention of diabetic renal disease with special reference to microalbuminuria. Lancet 346:1080-1084, 1995.
6. American Diabetes Association Diabetic nephropathy (Position Statement). Diabetes Care 20 (Suppl1.1):S24-S27, 1997
7. Hinberg, I. H, Katz, L, and Waddell, L, Sensitivity of in vitro diagnostic dipstick tests to urinary protein. Clin. Biochem, 11, 2, 63-64, 1978

SIMBOLI:

	Use by Date	IVD	In vitro diagnostic medical device
LOT	Batch Code	EC REP	Authorized Representative in the European Community
	Caution, consult Accompanying documents		Limit of Temperature
	Manufacturer	REF	Catalogue Number
	Contains sufficient for <n> tests	CONT	Contents
	Do not reuse	CE	Conformite Europeenne

Produttore:

 Shanghai Chemtron Biotech Co., Ltd.
NO.518, Qingdai Rd., International Medical Park,
Pudong District, Shanghai, P.R.China 201318
Tel: 86-21-68129666
Fax: 86-21-68129222
E-mail: international@chemtronbio.com
Web site: www.chemtronbio.com

EC REP European Representative:

Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestrasse 80, 20537 Hamburg Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

ITALY CONTACT

Meridian Healthcare srl
68, via G. Guglielmino
Tremestieri Etneo, CT 95030
Phone: +39 095 725 68 69
FAX: +39 095 725 44 54
www.meridianhealthcare.it