

G6PD QUALITATIVE TEST

Diagnostic Reagent kit for the *in vitro* determination of the activity of Glucose-6 Phosphate Dehydrogenase from the red cell hemolysate.



- Simple to perform dye reduction test.
- Available in convenient mono vial test pack.
- Visual interpretation of results, no need for any instruments.



ORDERING INFORMATION:

Product Code	Product	Pack-Size
96MB100	G6 PDH Visual test	10 test

G6PD TEST QUALITATIVO

Test qualitativo di screening a lettura visiva per la determinazione dell'attività della Glucosio 6 Fosfato Deidrogenasi (G6-PDH).

- Test di screening qualitativo.
- Reagente liquido e pronto all'uso.
- Metodo DCPIP 2,6-dichlorophenolindophenol .
- Monotest.
- Non richiede l'utilizzo di strumentazione.

