

HSV 1/2 ANTIGEN Rapid Test

StrongStep® HSV 1/2 antigen rapid test is a breakthrough advance in the diagnosis of HSV 1/2 for it is designated for qualitative detection of HSV antigen, which boasts of high sensitivity and specificity.

Benefits


- Simple, Accurate, Fast.
- Direct Testing on Mucocutaneous Lesions
- 🕒 Reliable Results in 10 minutes.
- Sensitivity : 89,10%
- Specificity: 95,90% Compared with PCR
- CE MARKED



HSV 1/2 ANTIGEN Rapid Test

StrongStep® HSV 1/2 è un test rapido lateral flow per la rilevazione dell'Antigene dell'HSV 1 e 2 da Lesioni Mucocutanee. HSV 1/2 Antigene è il primo test rapido alternativo al metodo PCR. Differenzia l'ulcera genitale da Herpes da altre tipologie di ulcera tipo la sifilide.

- Semplice, Accurato, Veloce.
- Test Diretto su tamponi vaginali, lesioni mucocutanee (ulcere, vescicole)
- 🕒 Risultato in 10 minuti.
- Sensibilità : 89,10%
- Specificità: 95,90% comparata con il metodo PCR
- CE MARKED



EUROPE Ordering Information:
HERPES ANTIGEN RAPID TEST
Product Code: 500070
Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)