


OnSITE® HBsAg Rapid Test

The *OnSite* HBsAg Combo Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of *hepatitis B surface antigen (HBsAg)* in human serum, plasma or whole blood.



- ✓ Simple, Accurate, Fast.
- ✓ Direct Testing in serum, plasma or whole blood.
 - 🕒 Reliable Results in 15 minutes.
 - ✓ Sensitivity : 100%
 - ✓ Specificity: 100%
- ✓ Analytic detection level as low as 0.5 ng/mL



ORDERING INFORMATION:



Code No.	Product	Pack-size
R0042C	OnSite® HBsAg Combo	30 Test

OnSITE® HBsAg Test Rapido

Test rapido qualitativo di screening per la determinazione del virus dell'Epatite B Antigene (HBsAg) in campioni di siero, plasma e sangue intero. Test ad esclusivo uso professionale in vitro.

- ✓ Semplice, Accurato, Veloce.
- ✓ Test diretto su siero, plasma o sangue intero.
- 🕒 Risultato in 15 minuti.
- ✓ Sensibilità: 100%.
- ✓ Specificità: 100%.
- ✓ Limite di rilevazione: 0.5 ng/mL

