

# OnSITE® HIV 1/2 Rapid Test

The *OnSite*® HIV 1/2 Ab Plus Combo Rapid Test is a lateral flow immunoassay for the simultaneous detection of anti-HIV-1 (including subtype O) and anti-HIV-2 antibodies (IgG, IgM, IgA) in human serum, plasma or whole blood. It is intended to be used by healthcare professionals as an aid in the diagnosis of infection with HIV. **3rd Generation, double-antigen based test to increase the sensitivity in early detection.**



For multiple screening



- ✓ Simple, Accurate, Fast.
- ✓ Direct Testing in serum, plasma or whole blood.
- 🕒 Reliable Results in 15 minutes.
- ✓ HIV-1 sensitivity : 100%
- ✓ HIV-2 sensitivity : 100%
- ✓ Specificity: 100%
- ✓ Performance validated with BBI low titre panel and FDA HIV Ab panel



#### ORDERING INFORMATION:

Code No.	Product	Pack-size
R0011C	<i>OnSite</i> ® HIV 1 & 2 Combo	30 Test

# OnSITE® HIV 1/2 Test Rapido

Test rapido di 3<sup>a</sup> generazione per la rilevazione differenziata degli anticorpi contro l'HIV (virus dell'immunodeficienza umana) di tipo 1 (compreso il sottotipo O) e di tipo 2 su campioni di sangue intero, siero o plasma.

- ✓ Semplice, Accurato, Veloce.
- ✓ Test diretto su siero, plasma o sangue intero.
- 🕒 Risultato in 15 minuti.
- ✓ Sensibilità: HIV-1 100%.
- ✓ Sensibilità: HIV-2 100%.
- ✓ Specificità: 100%.

