



25-OH-Vitamin D (total) ELISA

quantitative determination of total 25-OH Vitamin D (Vitamin D2 and Vitamin D3) in serum and plasma



25-OH-Vitamina D (total) ELISA

determinazione quantitativa della 25-OH Vitamina D totale (Vitamina D2 e Vitamina D3) su siero e plasma

25-OH Vitamin D ELISA KIT Vitamin D Total



SPECIFICHE

- ❖ Fase solida: enzyme-linked immunosorbent assay - ELISA
- ❖ Metodo: sequenziale competitivo
- ❖ Reagenti: liquidi e pronti all'uso
- ❖ Calibratori: 0, 5, 10, 25, 46, 85, 150 ng/ml
- ❖ Controlli: 1 alto e 1 basso inclusi nel kit
- ❖ Campione: siero
- ❖ Volume Campione: 10 µl
- ❖ Sensibilità: 1.2 ng/ml
- ❖ Effetto gancio su valori elevati: nessuno
- ❖ Interferenze da droghe o farmaci: nessuna
- ❖ Tempo di esecuzione saggio: 80 minuti
- ❖ Automatizzabile
- ❖ CE

25-OH Vitamin D ELISA KIT Vitamin D Total

FEATURES

- ❖ Name: 25-OH-Vitamin D (total) ELISA
- ❖ Description: ELISA for the quantitative determination of total 25-OH Vitamin D (Vitamin D2 and Vitamin D3) in serum
- ❖ Design: ELISA
- ❖ Sample Volume: 10 µl serum or plasma
- ❖ Total Assay Time: 80 min
- ❖ Calibrators: 0, 5, 10, 25, 46, 85, 150 ng/ml
- ❖ Controls: 2 vials
- ❖ Sensitivity: 1.2 ng/ml
- ❖ Storage: 2 - 8 °C
- ❖ Kit: 96 Determinations
- ❖ CE Mark: CE Marked



ORDERING INFORMATION:

25-OH-Vitamin D (total) ELISA

REF. 9425-300

96 TEST



related products:

Calcitonin U.S. ELISA

Calcitonin FAST ELISA

PTH Intact ELISA