

AMNIOQUICK® IGFBP-1 Rapid Test

Rapid immunochromatographic test, pen-shaped, for in vitro use only, dedicated to the detection of PROM (Premature Rupture Of foetal Membranes) by highlighting of the IGFBP-1 (Insulin Growth Factor Binding Protein-1)



- ◆ Simple, Accurate, Fast.
- ◆ Direct Testing from vaginal swab sample.
- ◆ Reliable Results in 10 minutes.
- ◆ Detection limit: 5 ng/mL.
 - ◆ Sensitivity : 100%
 - ◆ Specificity: 91%



ORDERING INFORMATION:

Code No.	Product	Pack-size
1090005	IGFBP-1 Amnioquick®	10 Test

AMNIOQUICK® IGFBP-1 Rapid Test

Test rapido per la determinazione della proteina IGFBP-1 (Insulin Growth Factor Binding Protein 1)
Rilevazione di eventuali rotture delle membrane fetali nella gestante.

- ◆ Semplice, Accurato, Veloce.
- ◆ Test Diretto da tampone vaginale.
- ◆ Risultato in 10 minuti.
- ◆ Positivo: Bande C e T presenza di liquido amniotico nel campione.
- ◆ Negativo: Banda C assenza di liquido amniotico nel campione.
- ◆ Sensibilità : 100%
- ◆ Specificità: 91%
- ◆ Limite di rilevazione: 5 ng/mL.

