





Smart Grip Research Micro Pipettes



A c c u r a c y a t m a x i m u m

Engineered to Perfection



1. **Soft Comfort**
 The smooth plunger mechanism offers least resistance for ease of effortless accuracy.
2. **Comfortable Finger Rest**
 for resting of the pipette on either hand.
3. **Tip ejector**
 for smooth ejection of the tip preventing any contamination.
4. **Ergonomics**
 Best in-hand grip for comfort.
5. **Precise setting**
 of the selected volume ensures accurate dispensing of the liquid.
6. **Light weight**
 of the pipette makes your job less tiresome and improves efficiency.
7. **Ideal Shaft Shape**
 Helps reach contents in narrow or deep tubes.
8. **Durable Tip Cone**
 offers highest resistance to shocks and chemical corrosion.

Calibration report enclosed with every pipette

Beautifully
 Color Coded









PRO Series Pipettes are having smart soft finger-grip of Thermoplastic Elastomer (TPE) ensuring

- Minimum transfer of body heat to the pipette resulting in high accuracy even on continuous use.
- Minimum fatigue of hands & fingers.
- Reduced thumb stress.
- Easy identification through color coding.
- Very attractive look.



Grip colour

					
VARIABLE 0.1 - 2.5 µl 0.5 - 10 µl 2 - 20 µl FIXED 2.5, 5, 10 µl MULTI 0.5 - 10 µl	VARIABLE 5 - 50 µl FIXED 20, 25 µl MULTI 5 - 50 µl	VARIABLE 10 - 100 µl FIXED 50 µl MULTI 10-100 µl	VARIABLE 20 - 200 µl 50 - 200 µl FIXED 100, 200 µl MULTI 20 - 200 µl	VARIABLE 100 - 1000 µl 200 - 1000 µl FIXED 250, 500 µl MULTI 30 - 300 µl	VARIABLE 0.5 - 5 ml, 1 - 5 ml 0.5 - 10 ml, 1 - 10 ml FIXED 1000, 2000 5000, 10000 µl



Smart Grip

Research Multichannel Pipettes



- Soft smart finger-grip of TPE.
- Manifold can rotate 360° for easy right or left hand operation.
- Easily removable & fully autoclavable manifold.
- One hand tip ejection.
- Soft force tip loading and ejection.
- Consistent sample loading.
- Leak free sealing of tip.
- Attractive color coding.
- Compatible with most types of tips.
- Calibration report enclosed with every pipette



wonderful accuracy & precision.

- Fully Autoclavable
- Color Coded
- Available in 5, 10, 20, 25, 50, and 100 µl fix volumes.

Multichannel Micropipettes

Cat. No.	Channel	Vol. Range µl
MAP 8 - 010	8 Ch.	0.5 - 10
MAP 8 - 050	8 Ch.	5 - 50
MAP 8 - 100	8 Ch.	10 - 100
MAP 8 - 200	8 Ch.	20 - 200
MAP 8 - 300	8 Ch.	30 - 300
MAP 12 - 010	12 Ch.	0.5 - 10
MAP 12 - 050	12 Ch.	5 - 50
MAP 12 - 100	12 Ch.	10 - 100
MAP 12 - 200	12 Ch.	20 - 200
MAP 12 - 300	12 Ch.	30 - 300

Micropipette Stand

Cat. No.	Capacity
SAP - 100	4 Micropipettes
SAP - 200	6 Micropipettes

Variable Volume Micropipettes

Cat. No.	Volume Range µl	Increment µl
VAP - 125	0.1 - 2.5	0.01
VAP - 100	0.5 - 10	0.1
VAP - 200	5 - 50	0.5
VAP - 300	50 - 200	1.0
VAP - 400	200 - 1000	5.0
VAP - 500	10 - 100	0.5
VAP - 600	100 - 1000	5.0
VAP - 700	2 - 20	0.1
VAP - 800	20 - 200	1.0
VAP - 900	1000 - 5000	50.0
VAP - 905	500 - 5000	50.0
VAP - 1000	1000 - 10000	100.0
VAP - 1005	500 - 10000	100.0

Fix Volume Micropipettes

Cat. No.	Volume Range µl
FAP - 0205	2.5
FAP - 0005	5
FAP - 0010	10
FAP - 0020	20
FAP - 0025	25
FAP - 0050	50
FAP - 0100	100
FAP - 0200	200
FAP - 0250	250
FAP - 0500	500
FAP - 1000	1000
FAP - 2000	2000
FAP - 5000	5000
FAP - 10000	10000

Mini Micropipettes

Meridian Healthcare srl

Via Caronda, 446 SC/A - 95129 Catania - Italy
Tel. +39 095 725 68 69 - Fax: +39 095 725 44 54
info@meridianhealthcare.it
www.meridianhealthcare.it



Fine Care Biosystems
www.accumaximum.com



ISO 17025 : 2005