



High Performance Digital Micro Volume Pipette

Meridian Healthcare[®]
medical diagnostics and equipment



Accuracy at maximum

Engineered to perfection

Soft Comfort

The smooth plunger mechanism offers least resistance for ease of operation.

Comfortable Finger Rest

for resting of the pipette on either hand.

Tip ejector

for smooth ejection of the tip preventing any contamination.

Precise Setting

of the selected volume with click-stop mechanism ensures accurate dispensing of the liquid.

Light weight

of the pipette makes your job less tiresome and improves efficiency.

Ideal Shaft Shape

helps reach contents in narrow or deep tubes.

Durable Tip Cone

offers highest resistance to shocks and chemical corrosion



- ◆ Calibration as per DIN 12650 & EN- ISO 8655 standards.
- ◆ " CE " mark awarded.
- ◆ Click - stop digital system for volume setting.
- ◆ Autoclavable lower assembly.
- ◆ Ergonomic Design ensuring light weight & soft plunger movement.
- ◆ High accuracy and precision guaranteed.
- ◆ Two step plunger operation allows "reverse pipetting technique".
- ◆ Individual serial No. permanently marked on each pipette.
- ◆ Calibration report enclosed with every pipette.
- ◆ Colour coded for easy identification.



Accumax

Multichannel Pipettes

- ◆ Soft force / leak free tip loading and ejection.
- ◆ Easily removable & fully autoclavable manifold.
- ◆ Easy Right or left hand operation because of 360° rotatable Manifold.
- ◆ One hand tip ejection.
- ◆ Attractive color coding.
- ◆ Compatible with most types of tips.



100% accuracy - 50% size - 100% efficiency

Beautifully Color Coded



VARIABLE	VARIABLE	VARIABLE	VARIABLE	VARIABLE	VARIABLE
0.1 - 2.5 μ l 0.5 - 10 μ l 2 - 20 μ l	5 - 50 μ l	10 - 100 μ l	20 - 200 μ l 50 - 200 μ l	100 - 1000 μ l 200 - 1000 μ l	0.5 - 5 ml 1 - 5 ml 0.5 - 10 ml 1 - 10 ml
FIXED	FIXED	FIXED	FIXED	FIXED	FIXED
2.5, 5, 10 μ l	20, 25 μ l	50 μ l	100, 200 μ l	250, 500 μ l	1, 2, 5, 10 ml
MULTI	MULTI	MULTI	MULTI	MULTI	
0.5 - 10 μ l	5 - 50 μ l	10 - 100 μ l	20 - 200 μ l	30 - 300 μ l	

Variable Vol. Micropipettes

Cat. No.	Vol. Range μ l	Increment μ l
VA - 125	0.1 - 2.5	0.01
VA - 100	0.5 - 10	0.1
VA - 200	5 - 50	0.5
VA - 300	50 - 200	1.0
VA - 400	200 - 1000	5.0
VA - 500	10 - 100	0.5
VA - 600	100 - 1000	5.0
VA - 700	2 - 20	0.1
VA - 800	20 - 200	1.0
VA - 900	1000 - 5000	50.0
VA - 905	500 - 5000	50.0
VA - 1000	1000-10000	100.0
VA - 1005	500-10000	100.0

Fix Vol. Micropipettes

Cat. No.	Vol. Range μ l
FA - 0205	2.5
FA - 0005	5
FA - 0010	10
FA - 0020	20
FA - 0025	25
FA - 0050	50
FA - 0100	100
FA - 0200	200
FA - 0250	250
FA - 0500	500
FA - 1000	1000
FA - 2000	2000
FA - 5000	5000
FA - 10000	10000

Multichannel Micropipettes

Cat. No.	Channel	Vol. Range μ l
MA 8 - 010	8 Ch.	0.5 - 10
MA 8 - 050	8 Ch.	5 - 50
MA 8 - 100	8 Ch.	10 - 100
MA 8 - 200	8 Ch.	20 - 200
MA 8 - 300	8 Ch.	30 - 300
MA 12 - 010	12 Ch.	0.5 - 10
MA 12 - 050	12 Ch.	5 - 50
MA 12 - 100	12 Ch.	10 - 100
MA 12 - 200	12 Ch.	20 - 200
MA 12 - 300	12 Ch.	30 - 300

Micropipette Stand

Cat. No.	Capacity
SA - 100	4 Micropipettes
SA - 200	6 Micropipettes

Meridian Healthcare srl

Tel. +39 095 725 68 69 Fax: +39 095 725 44 54
 info@meridianhealthcare.it
 www.meridianhealthcare.it

