

Trade diagnostical products  
for clinical laboratories  
products for life.

# Professional Product List 2024

## ALLERGY ELISA

The Specific IgE ELISA kit allows the quantitative detection of specific IgE in human serum.

A benefit of Specific IgE kit is the application of liquid biotinylated allergens binding monoclonal antibodies against human IgE immobilized on the inner surface of the wells. This technique has a number of advantages in comparison with the systems based on solid-phase allergens, particularly it avoids cross-reactions of anti-IgE antibodies with the immunoglobulins A, G, M and D.

## IgE SPECIFICHE

kit ELISA per la determinazione quantitativa delle IgE specifiche allergologiche nel siero umano.

Vantaggi del Kit IgE Allergen Basic: utilizzo di allergeni biotinilati liquidi leganti anticorpi monoclonali contro IgE umane, immobilizzati sulla superficie interna dei pozzetti. Nessuna cross-reaction di anticorpi anti-IgE con le immunoglobuline A, G, M e D.

Pagina lasciata  
intenzionalmente vuota

PRODUCT NAME	REF	PACKAGING
<b>ALLERGY - IgE SPECIFICHE</b>		
TOTAL IgE	2525-300	96 Test
AllergenP BASIC KIT	PA1000	96 Test
AllergenP BASIC KIT	PA1003	2x96 Test
AllergenP BASIC KIT	PA1004	5x96 Test
AllergenP BASIC KIT	PA1005	10x96 Test
AllergenP CALIBRATION SET	PA1100	10 curve
AllergenA SPECIFIC IgE POSITIVE CONTROL SERUM - INHALANTS	RA9249	10 Test
AllergenA SPECIFIC IgE POSITIVE CONTROL SERUM - FOODS	RA9250	10 Test
AllergenA SPECIFIC IgE NEGATIVE CONTROL SERUM	RA9251	10 Test
ALLERSCREEN INHALANTS MIX (D1, D2, E1, E2, E3, G2, G8, M3, M6, T4, T9, T11, W1, W6, W9, W21)	RAAX000	25 Test
INHALANTS MIX (D1, E1, E5, G6, G12, M2, T3, W6)	RAAX001	25 Test
ENVIRONMENT MIX (D1,D2,E1,E2)	RAAX002	25 Test
MITE MIX (D1, D2, D3, D70, D71, D72, D73, D74)	RADX002	25 Test
EPITHELIA MIX (E1, E2, E3, E4)	RAEX001	25 Test
EPITHELIA MIX (E1, E5, E6, E87, E88)	RAEX002	25 Test
RODENT MIX (E6, E82, E84, E87, E88)	RAEX070	25 Test
FEATHER MIX (E70, E85, E86, E89)	RAEX071	25 Test
DOMESTIC BIRDS FEATHERS MIX (E78, E93, E201, E213)	RAEX072	25 Test
EPITHELIA MIX (E1, E2, E3, E4, E5, E70, E81, E85, E86, E100)	RAEX100	25 Test
FOOD MIX (NUTS) (F13, F16, F17, F20, F36)	RAFX001	25 Test
FOOD MIX (SEA FOOD) (F3, F24, F37, F40, F41)	RAFX002	25 Test
FOOD MIX (CEREALS) (F4, F7, F8, F10, F11)	RAFX003	25 Test
FOOD MIX (SEA FOOD) (F3-F23-F24-F37)	RAFX004	25 Test
FOOD MIX (PAEDIATRIC) (F1, F2, F3, F4, F13, F14)	RAFX005	25 Test
FOOD MIX (SEA FOOD) (F3, F24, F37, F40, F41)	RAFX006	25 Test
FOOD MIX (CEREALS) (F4-F6-F7-F8-F9)	RAFX008	25 Test
FRUIT MIX (F20, F84, F87, F92, F259)	RAFX009	25 Test
FOOD MIX (PAEDIATRIC) (F1-F2-F3-F4-F14-F25-F75)	RAFX010	25 Test
VEGETABLE MIX (F12, F15, F31, F35)	RAFX013	25 Test
FOOD MIX (F25, F214, F216, F218)	RAFX014	25 Test
FOOD MIX (F33, F49, F92, F95)	RAFX015	25 Test
FOOD MIX (F44, F94, F208, F210)	RAFX016	25 Test
FRUIT MIX (F49, F92, F94, F95)	RAFX017	25 Test
FOOD MIX (CITRUS) F33, F208, F209, F302)	RAFX018	25 Test
FRUIT MIX (F84, F87, F92, F95, F210)	RAFX021	25 Test

FOOD MIX (CHEESE) (F70, F81, F82, F150, F198)	RAFX022	25 Test
SPICE MIX (F272, F273, F274, F275)	RAFX070	25 Test
SPICE MIX (F265, F267, F268, F282)	RAFX071	25 Test
SPICE MIX (F219, F269, F270, F271)	RAFX072	25 Test
FOOD MIX (FISH) (F3, F205, F206, F254)	RAFX074	25 Test
FOOD MIX (F1, F2, F4, F5, F8, F75, F76, F77, F78, F79, F81)	RAFX101	25 Test
FOOD MIX (F13, F14, F16, F17, F26, F45, F48, F83)	RAFX102	25 Test
FOOD MIX (F20, F25, F33, F44, F84, F87, F92, F95)	RAFX103	25 Test
STONE FRUITS MIX (F242, F95, F237, F255)	RAFX104	25 Test
FOOD MIX (F10, F12, F36, F84, F85, F93, F105, F221, F300)	RAFX105	25 Test
GRASS MIX (EARLY FLOWERING) (G2, G5, G6, G8, G10, G17)	RAGX001	25 Test
GRASS MIX (G3, G4, G5, G6, G8)	RAGX002	25 Test
GRASS MIX (LATE FLOWERING) (G1, G5, G6, G12, G13)	RAGX003	25 Test
DOG DANDER	RAE0005	25 Test
GUINEA PIG EPITHELIUM	RAE0006	25 Test
PIGEON DROPPINGS	RAE0007	25 Test
GOOSE FEATHERS	RAE0070	25 Test
MOUSE EPITHELIUM	RAE0071	25 Test
MOUSE URINE PROTEINS	RAE0072	25 Test
RAT URINE PROTEINS	RAE0074	25 Test
BUDGERIGAR DROPPING	RAE0077	25 Test
BUDGERIGAR FEATHERS	RAE0078	25 Test
BUDGERIGAR SERUM PROTEINS	RAE0079	25 Test
GOAT EPITHELIUM	RAE0080	25 Test
SHEEP EPITHELIUM	RAE0081	25 Test
RABBIT EPITHELIUM	RAE0082	25 Test
PIG EPITHELIUM	RAE0083	25 Test
HAMSTER EPITHELIUM	RAE0084	25 Test
CHICKEN FEATHERS	RAE0085	25 Test
DUCK FEATHERS	RAE0086	25 Test
RAT EPITH.+SERUM-URINE PROT.	RAE0087	25 Test
EPITH.+SERUM-URINE PROT. MOUSE	RAE0088	25 Test
TURKEY FEATHERS	RAE0089	25 Test
PARAKEET FEATHERS	RAE0093	25 Test
CAT DANDER	RAE0100	25 Test
CANARIAN FEATHERS	RAE0201	25 Test
NBOS D 6 BSA BOVINE SERUM ALBUMIN	RAE0204	25 Test
AFRICAN GREY PARROT FEATHERS	RAE0213	25 Test
FINCH FEATHERS	RAE0214	25 Test
PIGEON FEATHERS	RAE0215	25 Test
DERMATOPHAGOIDES PTERONYSSINUS	RAD0001	25 Test
DERMATOPHAGOIDES FARINAE	RAD0002	25 Test
DERMATOPHAGOIDES MICRO CERAS	RAD0003	25 Test
ACARUS SIRO	RAD0070	25 Test
LEPIDOGLYPHUS DESTRUCTOR	RAD0071	25 Test
TYROPHAGUS PUTRESENTIAE	RAD0072	25 Test
GLYCYPHAGUS DOMESTICUS	RAD0073	25 Test
EUROGLYPHUS MAYNEI	RAD0074	25 Test
BLOMIA TROPICALIS	RAD0201	25 Test
EGG WHITE	RAF0001	25 Test
COW MILK	RAF0002	25 Test

FISH (COD)	RAF0003	25 Test
WHEAT FLOUR	RAF0004	25 Test
RYE	RAF0005	25 Test
BARLEY	RAF0006	25 Test
OATS	RAF0007	25 Test
MAIZE	RAF0008	25 Test
RICE	RAF0009	25 Test
SESAME SEED	RAF0010	25 Test
BUCK WHEAT	RAF0011	25 Test
PEA	RAF0012	25 Test
PEANUT	RAF0013	25 Test
SOYA BEAN	RAF0014	25 Test
WHITE BEAN	RAF0015	25 Test
WALNUT	RAF0016	25 Test
HAZELNUT	RAF0017	25 Test
BRAZIL NUT	RAF0018	25 Test
ALMOND	RAF0020	25 Test
CRAB	RAF0023	25 Test
SHRIMP	RAF0024	25 Test
TOMATO	RAF0025	25 Test
PORK	RAF0026	25 Test
BEEF	RAF0027	25 Test
CARROT	RAF0031	25 Test
ORANGE	RAF0033	25 Test
TANGERINE	RAF0034	25 Test
POTATO	RAF0035	25 Test
COCONUT	RAF0036	25 Test
BLUE MUSSEL	RAF0037	25 Test
TUNA	RAF0040	25 Test
SALMON	RAF0041	25 Test
STRAWBERRY	RAF0044	25 Test
YEAST	RAF0045	25 Test
SWEET PEPPER	RAF0046	25 Test
GARLIC	RAF0047	25 Test
ONION	RAF0048	25 Test
APPLE	RAF0049	25 Test
MILLET	RAF0055	25 Test
DUCK MEAT	RAF0057	25 Test
OCTOPUS	RAF0059	25 Test
SARDINE	RAF0061	25 Test
KEFIR	RAF0063	25 Test
OVALBUMIN	RAF0067	25 Test
OVOMUCOID	RAF0068	25 Test
CONALBUMIN (OVOTRANSFERRIN)	RAF0069	25 Test
EGG YOLK	RAF0075	25 Test
ALPHA-LACTALBUMIN (BOVINE)	RAF0076	25 Test
BETA-LACTOGLOBULIN (BOVINE)	RAF0077	25 Test
CASEIN (BOVINE)	RAF0078	25 Test
GLUTEN	RAF0079	25 Test
LOBSTER	RAF0080	25 Test
CHEESE CHEDDAR TYPE	RAF0081	25 Test
CHEESE MOULD TYPE	RAF0082	25 Test
CHICKEN MEAT	RAF0083	25 Test
KIWI	RAF0084	25 Test
CELERY	RAF0085	25 Test
PARSLEY	RAF0086	25 Test
MELON	RAF0087	25 Test
MUTTON	RAF0088	25 Test
MUSTARD	RAF0089	25 Test
MALT	RAF0090	25 Test
MANGO	RAF0091	25 Test
BANANA	RAF0092	25 Test
COCOA	RAF0093	25 Test
PEAR	RAF0094	

PEACH	RAF0095	25 Test
AVOCADO	RAF0096	25 Test
CHOCOLATE	RAF0105	25 Test
RASPBERRY	RAF0111	25 Test
ZUCCHINI	RAF0113	25 Test
RED CABBAGE	RAF0116	25 Test
SPELT	RAF0124	25 Test
MINT	RAF0126	25 Test
PEPPER	RAF0125	25 Test
HIBISCUS	RAF0147	25 Test
CHEESE, EDAM TYPE	RAF0150	25 Test
ROQUEFORT CHEESE	RAF0152	25 Test
CHICORY	RAF0155	25 Test
RUCCOLA	RAF0171	25 Test
ARTICHOKE	RAF0172	25 Test
SEAWEED	RAF0173	25 Test
BREAM	RAF0185	25 Test
BRAN (WHEAT)	RAF0190	25 Test
QUAIL EGG	RAF0193	25 Test
CHEESE, GOUDA TYPE	RAF0198	25 Test
CASHEW NUT	RAF0202	25 Test
PISTACHIO NUT	RAF0203	25 Test
TROUT	RAF0204	25 Test
HERRING	RAF0205	25 Test
MACKEREL	RAF0206	25 Test
CLAM	RAF0207	25 Test
LEMON	RAF0208	25 Test
GRAPEFRUIT	RAF0209	25 Test
PINEAPPLE	RAF0210	25 Test
MUSHROOMS	RAF0212	25 Test
RABBIT MEAT	RAF0213	25 Test
SPINACH	RAF0214	25 Test
LETTUCE	RAF0215	25 Test
CABBAGE	RAF0216	25 Test
BRUSSEL SPROUTS	RAF0217	25 Test
SWEET PEPPER (SPICE)	RAF0218	25 Test
FENNEL SEED	RAF0219	25 Test
CIANNAMON	RAF0220	25 Test
COFFEE	RAF0221	25 Test
TEA	RAF0222	25 Test
POPPY (SEED)	RAF0224	25 Test
PUMPKIN	RAF0225	25 Test
RED RADISH	RAF0226	25 Test
CHINESE GREEN RADISH	RAF0227	25 Test
COOKED MILK	RAF0231	25 Test
GLIADIN	RAF0233	25 Test
LENTIL	RAF0235	25 Test
APRICOT	RAF0237	25 Test
CUTTLE FISH	RAF0239	25 Test
BEEF LIVER	RAF0241	25 Test
CHERRY	RAF0242	25 Test
CUCUMBER	RAF0244	25 Test
WHOLE EGG (HEN)	RAF0245	25 Test
GUAR GUM	RAF0246	25 Test
HONEY	RAF0247	25 Test
YOGURT	RAF0250	25 Test
CHEESE PARMESAN TYPE	RAF0251	25 Test
PINE-SEEDS	RAF0253	25 Test
PLAICE	RAF0254	25 Test
PLUM	RAF0255	25 Test
SQUID	RAF0258	25 Test
GRAPE	RAF0259	25 Test
BROCCOLI	RAF0260	25 Test
ASPARGUS	RAF0261	25 Test

EGG PLANT/AUBERGINE	RAF0262	25 Test
GREEN PEPPER	RAF0263	25 Test
EEL	RAF0264	25 Test
CARAWAY	RAF0265	25 Test
CARDAMON	RAF0267	25 Test
CLOVES	RAF0268	25 Test
BASIL	RAF0269	25 Test
GINGER	RAF0270	25 Test
ANISE	RAF0271	25 Test
TARRAGON	RAF0272	25 Test
THYME	RAF0273	25 Test
MAJORAN	RAF0274	25 Test
LOVAGE	RAF0275	25 Test
FENNEL	RAF0276	25 Test
DILL	RAF0277	25 Test
BAY LEAF	RAF0278	25 Test
CHILI PEPPER	RAF0279	25 Test
BLACK PEPPER	RAF0280	25 Test
CURRY	RAF0281	25 Test
NUTMEG	RAF0282	25 Test
TURKEY	RAF0284	25 Test
RED BEAN	RAF0287	25 Test
DATE	RAF0289	25 Test
OYSTER	RAF0290	25 Test
CAULIFLOWER	RAF0291	25 Test
GUAVA	RAF0292	25 Test
ARABIC GUM	RAF0297	25 Test
CHESTNUT	RAF0299	25 Test
GOAT MILK	RAF0300	25 Test
SPINY LOBSTER	RAF0304	25 Test
LIME	RAF0306	25 Test
HAKE	RAF0307	25 Test
SARDINE (PILCHARD)	RAF0308	25 Test
CHICK PEA	RAF0309	25 Test
SWORDFISH	RAF0312	25 Test
ANCHOVY	RAF0313	25 Test
SNAIL	RAF0314	25 Test
FRENCH BEAN	RAF0315	25 Test
CORIANDER	RAF0317	25 Test
CILANTRO	RAF0318	25 Test
BEET	RAF0319	25 Test
BLACKCURRANT	RAF0321	25 Test
RED CAVIAR	RAF0323	25 Test
HOPS	RAF0324	25 Test
EWE'S CHEESE	RAF0325	25 Test
GOAT CHEESE	RAF0326	25 Test
FIG	RAF0328	25 Test
WATERMELON	RAF0329	25 Test
SAFFRON	RAF0331	25 Test
CARP FISH	RAF0333	25 Test
SOLE	RAF0337	25 Test
OLIVE	RAF0342	25 Test
NECTARINE	RAF0343	25 Test
SAGE	RAF0344	25 Test
LYCHEE	RAF0348	25 Test
GROUPE	RAF0360	25 Test
HUNCHBACK SALMON	RAF0411	25 Test
ZANDER	RAF0415	25 Test
rApi g 1 CELERY	RAF0417	25 Test
ANTHOXANTHUM ODORATUM (SWEET VERNAL GRASS)	RAG0001	25 Test
CYNODON DACTYLON /(BERMUDA GRASS)	RAG0002	25 Test
DACTYLIS GLOMERATA (COCKSFOOT)	RAG0003	25 Test
FESTUCA ELATIOR (MEADOW FESCUE)	RAG0004	25 Test
LOLIUM PERENNE (RYE-GRASS)	RAG0005	25 Test

PHLEUM PRATENSE (TIMOTHY)	RAG0006	25 Test
PHRAGMITES COMMUNIS (COMMON REED)	RAG0007	25 Test
POA PRATENSIS (MEADOW GRASS)	RAG0008	25 Test
AGROSTIS STOLONIFERA (BENTGRASS)	RAG0009	25 Test
SORGHUM HALEPENSE (JOHNSON GRASS)	RAG0010	25 Test
BROMUS INERMIS (BROME GRASS)	RAG0011	25 Test
SECALE CEREALE (CULTIVATED RYE)	RAG0012	25 Test
HOLCUS LANATUS (VELVET GRASS)	RAG0013	25 Test
AVENA SATIVA (CULTIVATED OAT)	RAG0014	25 Test
TRITICUM AESTIVUM (CULTIVATED WHEAT)	RAG0015	25 Test
ALOPERCURUS PRATENSIS (MEADOW FOXTAIL)	RAG0016	25 Test
PASPALUM NOTATUM (BAHIA GRASS)	RAG0017	25 Test
AGROPYRON REPENS (COUCH (QUACK) GRASS)	RAG0021	25 Test
ELYMUS TRITICOIDES (WILD RYE GRASS)	RAG0070	25 Test
PHALARIS ARUNDINACEA (CANARY GRASS)	RAG0071	25 Test
HORDEUM VULGARE (BARLEY)	RAG0201	25 Test
ZEA MAYS (MAIZE)	RAG0202	25 Test
RPHL P1 PHLEUM PRATENSE (TIMOTHY)	RAG0205	25 Test
RPHL P2 PHLEUM PRATENSE (TIMOTHY)	RAG0206	25 Test
RPHL P12 PHLEUM PRATENSE (TIMOTHY)	RAG0212	25 Test
RPHL P 5 PHLEUM PRATENSE (TIMOTHY)	RAG0215	25 Test
PENICILLIUM NOTATUM	RAM0001	25 Test
CLADOSPORIUM HERBARUM	RAM0002	25 Test
ASPERGILLUS FUMIGATUS	RAM0003	25 Test
MUCOR RACEMOSUS	RAM0004	25 Test
CANDIDA ALBICANS	RAM0005	25 Test
ALTERNARIA ALTERNATA	RAM0006	25 Test
BOTRYTIS CINEREA	RAM0007	25 Test
HELMINTHOSPORIUM HALODES	RAM0008	25 Test
FUSARIUM MONILIFORME	RAM0009	25 Test
STEMPHYLIUM BOTRYOSUM	RAM0010	25 Test
RHIZOPUS NIGRICANS	RAM0011	25 Test
AUREOBASIDIUM	RAM0012	25 Test
TRICHODERMA VIRIDE	RAM0015	25 Test
CURVULARIA LUNATA	RAM0016	25 Test
ASPERGILLUS AMSTELODAMI	RAM0017	25 Test
PAECILOMYCES VARIOTII	RAM0024	25 Test
ASPERGILLUS NIGER	RAM0033	25 Test
SERPULA LACRYMANS	RAM0034	25 Test
SACCHAROMYCES CEREVISAE	RAM0044	25 Test
ASPERGILLUS ORYZAE	RAM0048	25 Test
TRICHOPHYTON RUBRUM	RAM0205	25 Test
MALASSEZIA SPP	RAM0227	25 Test
MOULD MIX (M1-M2-M3-M4-M6)	RAMX001	25 Test
JUTE	RAK0013	25 Test
SHEEP WOOL (TREATED)	RAK0020	25 Test
SILK (BOMBYX MORI)	RAK0074	25 Test
ISOCYANATE TDI	RAK0075	25 Test
ISOCYANATE MDI	RAK0076	25 Test
ISOCYANATE HDI	RAK0077	25 Test
ETHYLENE OXIDE	RAK0078	25 Test
PHTHALIC ANHYDRIDE	RAK0079	25 Test
FORMALDEHYDE	RAK0080	25 Test
FICUS BENJAMINA	RAK0081	25 Test
LATEX	RAK0082	25 Test
COTTON SEED	RAK0083	25 Test
SUNFLOWER SEED	RAK0084	25 Test
CHLORAMINE-T-HAS	RAK0085	25 Test
TRIMELLITIC ANHYDRIDE	RAK0086	25 Test
ALPHA-AMYLASE (FROM ASPERGILLUS ORYZAE)	RAK087A	25 Test
ALPHA-AMYLASE (FROM BARLEY MALT)	RAK087B	25 Test
HENNA	RAK0091	25 Test
PAPAIN-HAS	RAK0201	25 Test



BROMELAIN	RAK0202	25 Test
PHOSPHOLIPASE	RAK0203	25 Test
MEXATASE	RAK0204	25 Test
ALKALASE	RAK0205	25 Test
SAVINASE	RAK0206	25 Test
LYSOZYME	RAK0208	25 Test
PEPSIN	RAK0213	25 Test
BENZOIC ACID	RAK0300	25 Test
SORBIC ACID	RAK0301	25 Test
ASCARIS	RAP0001	25 Test
ANISAKIS	RAP0004	25 Test
ACER NEGUNDO (BOX-ELDER)	RAT0001	25 Test
ALNUS INCANA (GREY ALDER)	RAT0002	25 Test
BETULA ALBA (COMMON SILVER BIRCH)	RAT0003	25 Test
CORYLUS AVELLANA (HAZEL)	RAT0004	25 Test
FAGUS SPP (AMERICAN BEECH)	RAT0005	25 Test
JUNIPERUS SABINOIDES (JUNIPER)	RAT0006	25 Test
QUERCUS ROBUR (OAK)	RAT0007	25 Test
ULMUS SCABRA (ELM)	RAT0008	25 Test
OLEA EUROPEA (OLIVE)	RAT0009	25 Test
JUGLANS CALIFORNICA (WALNUT)	RAT0010	25 Test
PLATANUS ACERIFOLIA (LONDON PLANE, MAPLE LEAF SYCAMORE)	RAT0011	25 Test
SALIX CAPREA (WILLOW)	RAT0012	25 Test
POPULUS TREMULA (ASPEN)	RAT0013	25 Test
POPULUS DELTOIDES (COTTONWOOD)	RAT0014	25 Test
FRAXINUS EXCELSIOR (WHITE ASH)	RAT0015	25 Test
PINUS STROBUS (WHITE PINE)	RAT0016	25 Test
CRYPTOMERIA JAPONICA (JAPANESE CEDAR)	RAT0017	25 Test
EUCALYPTUS SPP (GUM-TREE)	RAT0018	25 Test
ACACIA SPP (ACACIA)	RAT0019	25 Test
PROSOPIS JULIFLORA (MESQUITE)	RAT0020	25 Test
MELALEUCA LEUCADENDRON (CAJEPUT TREE)	RAT0021	25 Test
CARYA PECAN (HICKORY)	RAT0022	25 Test
CUPRESSUS SEMPERVIRENS (ITALIAN CYPRESS)	RAT0023	25 Test
SYRINGA VULGARIS (LILAC)	RAT0024	25 Test
TILIA CORDATA (LINDEN)	RAT0027	25 Test
ROBINIA PSEUDOACACIA (FALSE ACACIA)	RAT0028	25 Test
PRUNUS AVIUM (SWEET CHERRY)	RAT0029	25 Test
MORUS ALBA (WHITE MULBERRY)	RAT0070	25 Test
CASUARINA EQUISETIFOLIA (AUSTRALIAN PINE)	RAT0073	25 Test
SALIX ALBA (WHITE WILLOW)	RAT0102	25 Test
PICEA EXCELSA (SPRUCE)	RAT0201	25 Test
AESCULUS HIPPOCASTANUM (HORSECHESTNUT)	RAT0203	25 Test
SAMBUCUS NIGRA (ELDER)	RAT0205	25 Test
CASTANEA SATIVA (CHESTNUT)	RAT0206	25 Test
PSEUDOTSUGA TAXIFOLIA (DOUGLAS FIR)	RAT0207	25 Test
CARPINUS BETULUS (EUROPEAN HORNBEAM)	RAT0209	25 Test
PHOENIX DACTYLIFERA (DATE PALME)	RAT0214	25 Test
RBET V1 BETULA VERRUCOSA (COMMON SILVER BIRCH)	RAT0215	25 Test
RBET V2 BETULA VERRUCOSA (COMMON SILVER BIRCH)	RAT0216	25 Test
RBET V4 BIRCH (BETULA VERRUCOSA)	RAT0220	25 Test
CUPRESSUS ARIZONICA (ARIZONA CYPRESS)	RAT0222	25 Test
AMBROSIA ELATIOR (COMMON RAGWEED)	RAW0001	25 Test
AMBROSIA PSILOSTACHYA (WESTERN RAGWEED)	RAW0002	25 Test
AMBROSIA TRIFIDA (GIANT RAGWEED)	RAW0003	25 Test
FRANSERIA ACANTHICARPA (FALSE RAGWEED)	RAW0004	25 Test
ARTEMISIA ABSINTHIUM (WORMWOOD)	RAW0005	25 Test
ARTEMISIA VULGARIS (MUGWORT)	RAW0006	25 Test
CHRISANTHEMUM LEUCANTHEMUM (OX-EYE DAISY)	RAW0007	25 Test
TARAXACUM VULGARE (DANDELION)	RAW0008	25 Test
PLANTAGO LANCEOLATA (PLANTAIN, RIBWORT)	RAW0009	25 Test
CHENOPODIUM ALBUM (GOOSEFOOT, LAMB'S QUARTERS)	RAW0010	25 Test
SALSOLA KALI (SALWORT)	RAW0011	25 Test
SOLIDAGO VIRGAUREA (GOLDEN ROD)	RAW0012	25 Test

XANTHIUM COMMUNE (COCKLEBUR)	RAW0013	25 Test
AMARANTHUS RETROFLEXUS (COMMON PIGWEED)	RAW0014	25 Test
ATRIPLEX LENTIFORMIS (SCALE)	RAW0015	25 Test
IVA CILIATA / IVA ANNUA (ROUGH MARSHELDER)	RAW0016	25 Test
KOCHIA SCOPARIA (FIREBUSH)	RAW0017	25 Test
RUMEX ACETOSELLA (SHEEP SORREL)	RAW0018	25 Test
PARIETARIA OFFICINALIS (WALL PELLITORY)	RAW0019	25 Test
URTICA DIOICA (NETTLE)	RAW0020	25 Test
PARIETARIA JUDAICA (WALL PELLITORY)	RAW0021	25 Test
PRIMULA VARIABILIS (PRIMEROSE)	RAW0036	25 Test
ALOE ARBORESCENS (KRANTS ALOE)	RAW0065	25 Test
BRASSICA NAPUS (RAPE)	RAW0203	25 Test
HELIANTHUS ANNUUS (SUNFLOWER)	RAW0204	25 Test
MATRICARIA CHAMOMILLA (CHAMOMILE)	RAW0206	25 Test
BETA VULGARIS (SUGAR BEET)	RAW0210	25 Test
NART V 1 ARTEMISIA VULGARIS (MUGWORT)	RAW0231	25 Test
NART V 3 ARTEMISIA VULGARIS (MUGWORT)	RAW0233	25 Test
PENICILLOYL G	RAC0001	25 Test
PENICILLOYL V	RAC0002	25 Test
ACTH	RAC0003	25 Test
CEFACLOR-HSA	RAC0007	25 Test
PARACETAMOL	RAC0020	25 Test
ASPIRIN-HAS	RAC0051	25 Test
PYRAZOLONE-HAS	RAC0052	25 Test
TRIMETHOPRIM-HAS	RAC0057	25 Test
SULFAMETHOXAZOLE-HSA	RAC0058	25 Test
ULTRACAINE (ARTICAINE)	RAC0059	25 Test
PHENYLBUTAZONE-HSA	RAC0065	25 Test
ULTRACAINE (ARTICAINE) -HSA	RAC0068	25 Test
INSULIN PORCIN	RAC0070	25 Test
INSULIN BOVINE	RAC0071	25 Test
INSULIN HUMAN	RAC0073	25 Test
GELATIN	RAC0074	25 Test
THEOPHYLLINE	RAC0081	25 Test
LIDOCAINE	RAC0082	25 Test
NOVOCAINE (PROCAINE)	RAC0083	25 Test
BENZOCAINE	RAC0084	25 Test
MEPIVACAINE	RAC0088	25 Test
PROPYPHENAZONE-HSA	RAC0090	25 Test
METAMIZOL-HSA	RAC0091	25 Test
INDOMETACIN-HSA	RAC0093	25 Test
NEOMYCIN	RAC0095	25 Test
4-AMINOBENZOIC ACID (VITAMIN B10)	RAC0105	25 Test
CAPTOPRIL	RAC0107	25 Test
CIPROFLOXACIN-HSA	RAC0108	25 Test
NAPROXEN	RAC0110	25 Test
OXACILLIN	RAC0116	25 Test
METRONIDAZOLE	RAC0153	25 Test
KETOPROFEN-HSA	RAC0172	25 Test
NORFLOXACIN-HSA	RAC0175	25 Test
EPINEPHRINE	RAC0196	25 Test
AMPICILLIN-HAS	RAC0203	25 Test
AMOXICILLIN-HAS	RAC0204	25 Test
TETRACYCLINE-HAS	RAC0205	25 Test
CEPHALOSPORIN-HAS	RAC0206	25 Test
GENTAMYCIN-HSA	RAC0207	25 Test
PROTAMINE	RAC0208	25 Test
CHYMOPAPAIN	RAC0209	25 Test
rPar j 2 PARIETARIA JUDAICA ( WALL PELLITORY)	RAC0211	25 Test
ERYTHROMYCIN-HAS	RAC0212	25 Test
DICLOFENAC	RAC0281	25 Test
IBUPROFEN	RAC0286	25 Test
STREPTOMYCIN-HSA	RAC0295	25 Test
RIFAMPICIN-HSA	RAC0301	25 Test

ACETYLCYSTEINE	RAC0320	25 Test
PREDNISOLONE-HSA	RAC0424	25 Test
SPYRAMICIN-HAS	RAC0436	25 Test
COTTON LINTERS	RAB0001	25 Test
BEECH WOOD	RAB0032	25 Test
OAK WOOD	RAB0033	25 Test
PINE WOOD	RAB0036	25 Test
ELM WOOD	RAB0049	25 Test
DAPHNIA	RAB0002	25 Test
TOBACCO DUST	RAB0003	25 Test
APIS MELLIFERA (HONEY BEE)	RAI0001	25 Test
DOLICHOVESPULA MACULATA (WHITE-FACED HORNET)	RAI0002	25 Test
VESPULA GERMANICA (COMMON WASP - YELLOW JACKET)	RAI0003	25 Test
POLISTES SPP (PAPER WASP)	RAI0004	25 Test
DOLICHOVESPULA ARENARIA (YELLOW HORNET)	RAI0005	25 Test
BLATELLA GERMANICA (COCKROACH)	RAI0006	25 Test
HETEROCERA MIX (MOTH)	RAI0008	25 Test
HOUSEFLY MUSCA DOMESTICA	RAI0015	25 Test
APHIDIDAE (APHIDS)	RAI0067	25 Test
SIMULIUM VENUSTRUM (BLACK FLY)	RAI0068	25 Test
FORMICA SPP.( RED WOOD ANT)	RAI0069	25 Test
SOLENOPSIS INVICTA (FIRE ANT)	RAI0070	25 Test
AEDES COMUNIS (COMMON MOSQUITO)	RAI0071	25 Test
CULEX PIPIENS (MOSQUITO)	RAI0074	25 Test
VESPA CRABRO (EUROPEAN HORNET)	RAI0075	25 Test
PERIPLANETA AMERICANA (AMERICAN COCKROACH)	RAI0206	25 Test

ANALIZZATORE ELISA **ALISEI** CON  
PROGRAMMA ALLERGIE



Utilizzo di allergeni biotinilati liquidi leganti anticorpi monoclonali contro IgE umane, immobilizzati sulla superficie interna dei pozzetti. Nessuna cross-reaction di anticorpi anti-IgE con le immunoglobuline A, G, M e D.


Pagina lasciata  
intenzionalmente vuota

# CERTIFICAZIONI

Tutti i test presenti in questo catalogo sono certificati CE e soddisfanno tutte le disposizioni della direttiva 98/79/CE e successive modifiche ed integrazioni che la riguardano - salvo diversamente indicato.

## CE marking of In Vitro Diagnostic Medical Devices - IVDD (98/79/EC)

According to the IVDD, in vitro diagnostic medical devices include: reagents, reagent products, calibration materials, control materials, kits, instruments, apparatus, equipment, and systems that are intended for use in the examination of specimens taken from the human body (tissue, blood, urine, etc.) to diagnose diseases, to monitor a person's state of health, or to monitor therapeutic procedures.



Always There

Meridian Healthcare®  
medical diagnostics and equipment



**Italy Corporate Headquarter**  
Meridian Healthcare  
Via G. Guglielmino, 68  
95030 Tremestieri Etneo (Catania) – Italy  
Tel: +39 095 7256869 / +39 0922 83 21 88  
Fax: +39 095 725 44 54

Technical Support  
Engineering Department  
Meridian Healthcare  
Tel: +39 393 84 219 67  
Fax: +39 095 725 44 54  
mail: specialist@meridianhealthcare.it

**Customer Service**  
Meridian Healthcare  
Tel: +39 095 725 68 69 / +39 0922 83 21 88  
Fax: +39 095 725 44 54  
mail: info@meridianhealthcare.it  
marketing@meridianhealthcare.it

**Sales and Orders**  
Meridian Healthcare  
Tel: +39 095 725 68 69  
Fax: +39 095 725 44 54  
mail: ordini@meridianhealthcare.it  
amministrazione@meridianhealthcare