

## BACTERIAL VAGINOSIS Rapid Test

The NADAL pH Test is for measuring the pH value of vaginal fluid. The vaginal pH value is a useful screening parameter for disturbances of the vaginal microflora. Abnormal pH values are frequently an indicator for clinical or sub-clinical vaginal infections.

- Easy to use: can be used a self-testing for patients .
- Pain-Free method of determining the vaginal pH level.
- Reliable, Hygienic.
- ⌚ Reliable Results in 1 minute.
- Sensitivity : 99.9%
- Specificity: 86.4%



### ORDERING INFORMATION:



Product Code	Product	Pack-Size
194002	BACTERIAL Vaginosis Rapid Test	50 test



## VAGINOSI BATTERICA Test Rapido

Il test NADAL pH è concepito per la misurazione dei livelli di pH vaginale. La misurazione del pH vaginale è uno strumento di screening utile per l'individuazione di disturbi della microflora vaginale. Valori anormali di pH sono generalmente indicativi della sussistenza di infezioni vaginali cliniche o subcliniche. Test diretto con applicatore e scala colorimetrica.

- Semplice utilizzo anche in self-test.
- Affidabile, Igienico e Indolore.
- Applicatore incluso
- ⌚ Risultato in 1 minuto.
- Sensibilità : 99.9%
- Specificità : 86.4%



Lettura

