

BENCE JONES dipstick Rapid Test

Competitive immunochromatographic rapid test Vertical Flow for detection of free light chains K and λ in urine.



- Simple, Accurate, Fast.
- Direct Testing on urine or urines 24 hours collection .
 - 🕒 Reliable Results in 15 minutes.
- Kit sensitivity for K and λ free light chains is 30 mg/L.
 - Specificity: 100%
 - Precision: 100%



ORDERING INFORMATION:

| Product Code | Product | Pack-Size |
|--------------|---------------------------------|-----------|
| MH81200 | BJ Free k/λ Dipstick Rapid Test | 25 test |

BENCE JONES dipstick Rapid Test

Test rapido immunocromatografico Vertical Flow a competizione per la determinazione delle catene leggere libere K e λ nelle urine.

- Semplice, Accurato, Veloce.
- Test diretto su urine o urine delle 24 h.
- 🕒 Risultato in 15 minuti.
- Sensibilità : per le catene leggere libere K e λ è di 30 mg/l.
- Specificità : 100%
- Precisione: 100%

