

Lyme Borreliosis



Lyme Borreliosis

Test Lateral Flow for detection of human IgG and IgM to European strains of *Borrelia* in human serum and plasma

LATERAL FLOW

DETECTION

Specific Antibody IgG & IgM


Rapid Test Visual

Sensitivity: 91,3%

Specificity: 81,3%

Detection Time: 10 minutes

CE marked



CE

ORDERING INFORMATION:

REF.
870001

Product
Borrelia Burgdorferi Rapid Test

Packaging
10 test