

# CLOSTRIDIUM DIFFICILE Rapid Test

*C. difficile* GDH + Toxin A + Toxin B is a rapid chromatographic immunoassay combo card for the simultaneous qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A and Toxin B in human faeces that aids in the diagnosis of *C.difficile* infection. AVAILABLE: **TOXIN A + B** - only **GDH** and **A + B + GDH** combo.

## FEATURES

- 🕒 Simple, Accurate, Fast.
- Direct Testing on human feces specimens.
- Ambulatory can be used.
- 🕒 Reliable Results in 10 minutes.
- Sensitivity A+B: 98,2%
- Specificity A+B: 98,6%
- Sensitivity GDH: >99,0%
- Specificity GDH: >99,0%
- CE marked



JusChek®

# CLOSTRIDIUM DIFFICILE Rapid Test

*C. difficile* GDH + Toxin A + Toxin B è un test rapido immunocromatografico per la rilevazione simultanea e differenziata di *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A and Toxin B su campioni fecali come ausilio per la diagnosi di infezioni da *C.difficile*. DISPONIBILI: **TOXIN A + B** - solo **GDH** e **A + B + GDH** combo.

## SPECIFICHE

- 🕒 Semplice, Accurato, Veloce.
- Test diretto su campioni di feci.
- Utilizzo previsto anche in ambulatorio.
- 🕒 Risultato in 10 minuti.
- Sensibilità Tossine A+B: 98,2%
- Specificità Tossine A+B: 98,6%
- Sensibilità GDH: >99,0%
- Specificità GDH: >99,0%
- CE marked



## EUROPE Ordering Information:



### PRODUCT

CLOSTRIDIUM Difficile TOXINA A + B

CLOSTRIDIUM Difficile Ag (GDH)

CLOSTRIDIUM Difficile TOXINA A + B + GDH

### REF.

ICDT-625

ICDG-602

ICDT-635

### N. Test

10 card

10 card

10 card

