

## **CALCITONIN FAST ELISA**

**Calcitonin FAST: For the Diagnosis of Medullary Carcinoma of the Thyroid and Neoplastic Disorders.**

Meridian Healthcare®

**elisa calcitonin FAST**

## **CALCITONINA FAST ELISA**

**determinazione quantitativa della Calcitonina su siero. Metodo FAST 105 minuti.**

## CALCITONINA Fast

### ELISA KIT Calcitonina Fast



## CALCITONIN Fast

### ELISA KIT Calcitonin Fast

#### SPECIFICHE


- ❖ Fase solida: enzyme-linked immunosorbent assay - ELISA
- ❖ Metodo: Sandwich Equilibrium Method
- ❖ Reazione: colorimetrica direttamente proporzionale
- ❖ Controllo: incluso nel kit
- ❖ Calibratori: 0, 10, 40, 150, 400, 1000 pg/ml
- ❖ Campione: siero
- ❖ Lettura: 450 nm o 405 nm per conc >300 pg/ml.
- ❖ Sensibilità: 1 pg/ml
- ❖ Tempo di esecuzione saggio: 105 minuti
- ❖ CE Marked

#### FEATURES

- ❖ Name: Calcitonin Fast ELISA
- ❖ Description: Microplate Enzyme Immunoassay, Colorimetric
- ❖ Design: ELISA
- ❖ Sample Volume: 50 µl serum
- ❖ Total Assay Time: 105 minutes
- ❖ Sensitivity: 1 pg/ml.
- ❖ Calibrators: 0, 10, 40, 150, 400, 1000 pg/ml Dried
- ❖ Reading: 450 nm or 405 nm for conc >300 pg/ml
- ❖ Kit: 96 Determinations
- ❖ CE Mark: CE Marked

# CALCITONIN

## FAST ELISA



**ORDERING INFORMATION:**  
**CALCITONIN FAST ELISA**  
**REF. 9325-300A**  
**96 TEST** 

**related products:**

**25-OH-Vitamin D (total) ELISA** 

**PTH Intact ELISA** 