




CALCITONIN FAST ELISA

Calcitonin FAST: For the Diagnosis of Medullary Carcinoma of the Thyroid and Neoplastic Disorders.

Meridian Healthcare®

elisa calcitonin FAST

CALCITONINA FAST ELISA

determinazione quantitativa della Calcitonina su siero. Metodo FAST 105 minuti.

CALCITONINA Fast

ELISA KIT Calcitonina Fast



CALCITONIN Fast

ELISA KIT Calcitonin Fast

SPECIFICHE


- ❖ Fase solida: enzyme-linked immunosorbent assay - ELISA
- ❖ Metodo: Sandwich Equilibrium Method
- ❖ Reazione: colorimetrica direttamente proporzionale
- ❖ Controllo: incluso nel kit
- ❖ Calibratori: 0, 10, 40, 150, 400, 1000 pg/ml
- ❖ Campione: siero
- ❖ Lettura: 450 nm o 405 nm per conc >300 pg/ml.
- ❖ Sensibilità: 1 pg/ml
- ❖ Tempo di esecuzione saggio: 105 minuti
- ❖ CE Marked

FEATURES

- ❖ Name: Calcitonin Fast ELISA
- ❖ Description: Microplate Enzyme Immunoassay, Colorimetric
- ❖ Design: ELISA
- ❖ Sample Volume: 50 µl serum
- ❖ Total Assay Time: 105 minutes
- ❖ Sensitivity: 1 pg/ml.
- ❖ Calibrators: 0, 10, 40, 150, 400, 1000 pg/ml Dried
- ❖ Reading: 450 nm or 405 nm for conc >300 pg/ml
- ❖ Kit: 96 Determinations
- ❖ CE Mark: CE Marked

CALCITONIN

FAST ELISA



ORDERING INFORMATION:
CALCITONIN FAST ELISA
REF. 9325-300A
96 TEST



related products:

25-OH-Vitamin D (total) ELISA



PTH Intact ELISA

