

# Creatine Kinase MB (CK-MB) Rapid Test

The Meridian Healthcare CK-MB Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of Creatine Kinase MB in human whole blood, serum, and plasma samples. This test is a tool to assist in the diagnosis of myocardial infarction.


- Simple, Accurate, Fast.
- Sample: serum, plasma and W.B.
- ⌚ Reliable Results in 10 minutes.
- Sensitivity : 99.9%
- Specificity: 99.4%
- Cut-Off: 5 ng/ml.



# Creatina Chinasi MB (CK-MB) Rapid Test

Test rapido immunocromatografico a flusso laterale per la determinazione qualitativa della Creatinchinasi MB (CK-MB) nel siero, plasma e sangue intero.

- Semplice, Accurato, Veloce.
- Test su siero, plasma e sangue intero.
- ⌚ Risultato in 10 minuti.
- Sensibilità : 99.9%
- Specificità : 99.4%
- Cut-Off: 5 ng/ml.



Positive

Negative

Invalid

EUROPE Ordering Information:  
**CK - MB RAPID TEST**  
Product Code: **RAPG-CKMB-001**  
Packaging: **10 test**

Meridian Healthcare®

Meridian Healthcare srl  
68 - Tremestieri Etneo  
Catania, IT 95030  
Phone: 095 725 6869 (0039)  
Fax: 095 725 4454 (0039)