

# CRYPTOCOCCUS Antigen Rapid Test

The product is an immunochromatographic test system for the qualitative detection of the capsular polysaccharide antigens of Cryptococcus species complex in human serum, and cerebrospinal fluid (CSF). It is a single-use laboratory test and can be used as an aid in the diagnosis of cryptococcosis.

## Features

- Simple, Accurate, Fast
- Direct Testing on Serum and CSF
- ⌚ Reliable Results in 10 minutes
  - Sensitivity : 99.60%
  - Specificity: 99.78%
- Cut-off value: 0.5 ng/mL.
  - CE marked



# CRYPTOCOCCUS Antigen Test Rapido

Nuovo dispositivo lateral flow per la rilevazione di antigeni polisaccaridi associati a un'infezione derivante da Cryptococcus neoformans e Cryptococcus gattii. Test qualitativo a cut-off, di facile esecuzione e della durata di 10 minuti. Come campioni possono essere utilizzati siero o liquido cerebrospinale (CSF).

## Specifiche

- Semplice, Accurato, Veloce
- Test diretto su siero o CSF
- ⌚ Risultato in 10 minuti
- Sensibilità : 99.60%
- Specificità: 99.78%
- Cut-off: 0.5 ng/mL.
- CE marked



### Controlli disponibili:

- Ctrl Positivo: ordinabile separatamente

**EUROPE Ordering Information:**  
**CRYPTOCOCCUS ANTIGEN RAPID TEST**  
Product Code: BIO-1125  
Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl  
68 - Tremestieri Etneo  
Catania, IT 95030  
Phone: 095 725 6869 (0039)  
Fax: 0922 83 2188 (0039)

