

D-DIMER Lateral Flow

The D-dimer Rapid Test Cassette (Whole Blood/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of human D-dimer in whole blood or plasma as an aid in the diagnosis of Disseminated Intravascular Coagulopathy (DIC), deep venous thrombosis (DVT) and pulmonary embolism (PE).

FEATURES:

- Simple, Accurate, Fast.
- Sample: whole blood or plasma.
- 🕒 Reliable Results in 10 minutes.
- Limit of detection: 500 ng/mL.
- Sensitivity: 96.5%
- Specificity: 98.3%



D-DIMERO Lateral Flow

Dispositivo medico-diagnostico in vitro, su supporto immunocromatografico rapido per la rilevazione qualitativa del D-dimero su sangue intero (prelievo venoso e da puntura del polpastrello), siero o plasma da utilizzare come ausilio nella diagnosi di trombosi venosa profonda (TVP), embolia polmonare (EP) e coagulazione intravascolare disseminata (CID). Esclusivamente per uso diagnostico professionale in vitro.

SPECIFICHE:

- Campione: sangue intero o plasma.
- 🕒 Risultato in 10 minuti.
- Limite di rilevazione: 500 ng/mL.
- Sensibilità: 96.5%
- Specificità: 98.3%



EUROPE Ordering Information:

D-DIMER RAPID TEST
Product Code: CDM-402
Packaging: 10 test

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)