

DIP-SLIDE

Dip-Slide for the microbiological screening in **Urinary Tract Infections**



Description	Packaging	Ref.
URITEST (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500152 51015
URITEST N (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides 500 slides	500232 51023 51123
URITEST C (MAC CONKEY MUG AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500242 51024
URITEST EC (MAC CONKEY AGAR / E.COLI AGAR / CLED AGAR) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20 slides 120 slides	500412 51041
URITEST 2 (CLED AGAR/ MAC CONKEY AGAR) Total bacterial count, enterobacteria.	20 slides 120 slides 500 slides	500302 51030 51130
URITEST M (MAC CONKEY AGAR / MALT SELECTIVE AGAR / CLED ANDRADE AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500182 51018 51118
URITEST MALTO (MAC CONKEY AGAR / MALT AGAR / CLED AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500402 51040 51140
URITEST EF (MAC CONKEY AGAR / BILE ESCULIN AGAR / CLED AGAR) Enterobacteria, enterococci, total bacterial count.	20 slides 120 slides	500702 51070
URITEST PENTA (CLED AGAR / MAC CONKEY AGAR / CETRIMIDE AGAR / MALT AGAR / BILE ESCULIN AGAR) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20 slides 120 slides	500142 51014

Dip-Slide for the microbiological screening in **Urogenital Infections**

Description	Packaging	Ref.
VAGITEST (CHOCOLATE / ROGOSA / THAYER MARTIN) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10 slides 120 slides	50020 51020

Dip-Slide for the microbiological screening in **Skin Infections**

Description	Packaging	Ref.
DERMATEST (MALTO / CORN MEAL / DTM-TAPLIN) Pathogenic dermatophytes and others non pathogenic fungi.	10 slides 120 slides	50021 51021

ITALY CONTACT:

Meridian Healthcare srl

Via G. Guglielmino, 68 - 95030 Tremestieri Etneo -

CT Tel. 095 725 68 69 - Fax: 095 725 44 54

e-mail: info@meridianhealthcare.it

Meridian Healthcare®