

EmA ENDOMYSIUM IgA ANTIBODIES

The IFA test EmA IgA is a immunofluorescence Assay for the semiquantitative detection of IgA specific anti-endomysial antibodies in human serum and plasma.



- Standardized test format.
- Ready to use reagents.
- Teflon mask slides ensure easy handling during test procedure.
- Individually sealed and vacuum packed single slides for a long shelf-life.
- Positive Control: ready to use
- Negative Control: ready to use
 - Sensibility: 100%
 - Specificity: 99,3%
 - Test: 48 or 96

ORDERING INFORMATION:



Code No.	Product	Pack-size
EDM4K	Anti-Endomysium IgA	12x4 wells
EDM8K	Anti-Endomysium IgA	12x8 wells

EmA anti-ENDOMISIO IgA IFA

Test in immunofluorescenza per la rilevazione semiquantitativa degli anticorpi anti-endomisio IgA specifici nel siero e plasma umano.

- Formato del test standardizzato.
- Reagenti liquidi e pronti all'uso.
- Slides teflonati.
- Controllo Positivo: liquido e pronto all'uso.
- Controllo Negativo: liquido e pronto all'uso.
- N. Test: 48 o 96
- Sensibilità: 100%
- Specificità: 99,3%

