

ENTEROVIRUS Ag VP1 PEPTIDE Rapid Test

The Rapid-VIDITEST Enterovirus Card is a rapid chromatographic immunoassay for the qualitative detection of Enterovirus antigens (VP1 peptide) in faecal samples to aid in the diagnosis of Enterovirus infection. Only for laboratory use.



Features

- Simple, Accurate, Fast
- Direct Testing on Faeces
- Ambulatory can be used
- 🕒 Reliable Results in 10 minutes
 - Sensitivity : >99%
 - Specificity: >99%
 - CE marked



ENTEROVIRUS Ag VP1 PEPTIDE Test Rapido

Metodo sandwich immunocromatografico che impiega anticorpi monoclonali coniugati con latex colorato, specifici per gli antigeni di Enterovirus e anticorpi monoclonali su fase solida, per la determinazione diretta di Enterovirus (Poliovirus, Coxsackievirus, Echovirus: sono virus RNA a singolo filamento) nelle feci, con elevati livelli di sensibilità e specificità.

Specifiche

- Semplice, Accurato, Veloce
- Test diretto su campioni fecali
- Utilizzo in laboratorio e in ambulatorio
- 🕒 Risultato in 10 minuti
- Sensibilità : >99%
- Specificità: >99%
- CE marked

Controlli disponibili:

- Ctrl Positivo: ODZ-203 kit da 10 test



EUROPE Ordering Information:
ENTEROVIRUS Ag VP1 Peptide RAPID TEST
Product Code: ODZ-187
Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)

