

TEST GLUTEN Deamidated Gliadin IgG Rapid Test

EXACTO® Test Gluten is an easy-to-use, reliable self-test that can be performed at home for the initial screening of celiac disease, also known as gluten intolerance. **EXACTO® Test Gluten** is an immunochromatographic test for the detection of antibodies specifically produced in case of celiac disease and directed against deamidated gliadin peptides.



PERFORMANCES

- Simple, Accurate, Fast.
- Direct Testing in whole blood.
- 🕒 Reliable Results in 10 minutes.
 - Sensitivity : 93.8%
 - Specificity : 97.9%

KIT Content

- 1 Test
- 1 Alcohol wipe
- 1 Lancet
- 1 Plaster



Gliadina Deamidata IgG Rapid Test

EXACTO® Test Gluten è un test immunocromatografico Lateral Flow per la rilevazione di anticorpi specificamente prodotti in caso di celiachia e diretti contro i peptidi di gliadina deamidata.

SPECIFICHE

- ✓ Semplice, Accurato, Veloce.
- ✓ Test diretto su sangue intero.
- 🕒 Risultato in 10 minuti.
- ✓ Sensibilità : 93.8%
- ✓ Specificità : 97.9%
- ✓ Kit completo di lancette pungidito.



Codice Dispositivo
858865

Dispositivo
EXACTO® Test Gluten IgG Rapid Test

Confezione
1 test