

Meridian Healthcare®

D I A G N O S T I C S

tHb


K⁺

pH



CO₂

Glu



Lac

Na



MISSION®

Quality Controls for Blood Gas, ISE, Coox and Metabolite systems



MISSION®



Mission CONTROL™

Aqueous Based Controls

- Assayed for pH, pCO₂, pO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, Li⁺, tCO₂
- 36 month stability room temperature
- 60 second Open Vial Stability

Assayed Instrumentation

Abbott I-Stat
 Idexx
 Instrumentation Laboratory
 ITC (formerly Diametrics)
 Medica
 Nova Biomedical
 OPTI Medical
 Radiometer
 Roche Diagnostics (formerly AVL)
 Siemens Diagnostics (formerly Bayer)

Mission TRINITY™

Aqueous Based Control with Dyes and Metabolites

- Assayed for pH, pCO₂, pO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, Li⁺, Glucose, Lactate, tHb, sO₂, MetHb, COHb, HHb, O₂Hb
- 36 month stability when refrigerated
- 60 second Open Vial Stability

Assayed Instrumentation

Instrumentation Laboratory
 Radiometer
 Roche Diagnostics (formerly AVL)
 Siemens Diagnostics (formerly Bayer)

Mission COMPLETE™

Aqueous Based CVM Linearity Control

- Assayed for pH, pCO₂, pO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, Li⁺, Glucose, Lactate
- 36 month stability when refrigerated
- 60 second Open Vial Stability
- 5 Levels with Color Code Ampules covering analyzer significant ranges

Assayed Instrumentation

Abbott I-Stat
 Instrumentation Laboratory
 ITC (formerly Diametrics)
 Nova Biomedical
 OPTI Medical
 Radiometer
 Roche Diagnostics (formerly AVL)
 Siemens Diagnostics (formerly Bayer)



Meridian Healthcare®

Meridian Healthcare srl

Tel. +39 095 725 68 69 Fax: +39 095 725 44 54
 info@meridianhealthcare.it
 www.meridianhealthcare.it

Meridian Healthcare®