

## FTA-ABS (Syphilis) Titratable IFA Kit

The Diagnostic Automation, Inc. FTA-ABS IFA test system is designed for the qualitative determination of antibodies of *Treponema pallidum* (T. Pallidum), and is intended to be used as an aid in the confirmation of syphilis antibodies.



- FTA-Abs Titratable IFA Kit
- Test: 10 x 10 wells
- ⌚ Reliable Results in 110 minutes.
- Principle: Quantitative
- Sample: 50 µL serum
- Method: Fluorescent Treponemal Antibody Absorption
- Shelf Life: 12 Months

### ORDERING INFORMATION:



Product Code  
351010-T

Product  
FTA-ABS (Syphilis) IFA

Pack-Size  
10x10 test

## FTA-ABS (Syphilis) IFA Kit Quantitativo

Fluorescent Treponemal Antibody-Absorption test immunologico in immunofluorescenza specifico utilizzato per la diagnosi di infezione luetica, o sifilitica. Il dispositivo rileva la presenza di anticorpi diretti contro il *Treponema pallidum*, agente causale della sifilide.

- Test IFA Quantitativo
- Numero Test: 10 x 10 pozzetti
- Campione: Siero - 50 µL.
- ⌚ Risultato in 110 minuti
- Emivita: 12 mesi

### Produttore

Diagnostic Automation  
Cortez Diagnostics, Inc. 21250 Califa Street, Suite 102 and 116,  
Woodland Hills, California 91367 USA

