

OnSite™ EPATITE A IgG/IgM

OnSite HAV IgG/IgM Rapid Test is the most sensitive HAV rapid test on the market for quick and accurate determination of Hepatitis A infection in a single test at every stage: early acute infection (IgM positive) - late acute infection (IgG and IgM positive) - protective immunity (IgG positive).

FEATURES

- Simple, Accurate, Fast.
- Sample: serum, plasma or W.B.
- ⌚ Reliable Results in 15 minutes.
- Sensitivity IgG: 100%
- Specificity IgG: 94.7%
- Sensitivity IgG: 94.8%
- Specificity IgG: 96.7%



OnSite™ EPATITE A IgG/IgM

Nuovo Test Rapido Lateral Flow a doppia banda differenziata per lo screening rapido qualitativo degli anticorpi anti epatite A IgG e IgM in campioni di siero, plasma e sangue intero. Test ad esclusivo uso professionale in vitro.

SPECIFICHE

- ◆ Test a doppia banda differenziata per IgG e IgM
- ◆ Campione: siero, plasma e sangue intero
- ⌚ Tempo del saggio: 15 minuti
- ◆ Sensibilità IgG: 100% - Specificità IgG: 94.7%
- ◆ Sensibilità IgM: 94.8% - Specificità IgM: 96.7%



EUROPE Ordering Information:
HAV IgG/IgM RAPID TEST
Product Code: R0092C
Packaging: 30 test

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)