

HEPATITIS DISEASE

HBsAg

IMMUNOASSAY LATERAL FLOW

AllTest™ HBsAg

HBsAg Lateral Flow

Test rapido qualitativo di screening per la determinazione del virus dell'Epatite B Antigene (HBsAg) in campioni di siero, plasma e sangue intero. Test ad esclusivo uso professionale in vitro.

SPECIFICHE:

- ✓ Test diretto su siero e plasma e sangue intero e capillare.
- 🕒 Tempo di esecuzione: 15 min.
- ✓ Campione: 3 gocce siero o plasma / 3 gocce di sangue intero o da polpastrello.
- ✓ Sensibilità: 99.9%
- ✓ Specificità: 99.9%

Lancette sterili e Alcohol pad possono essere forniti separatamente.

INFO PER ORDINI:

| REF | Product | Pack-size |
|-----------|------------------|-----------|
| IHBSG-402 | HBsAg Rapid Test | 25 test |



AllTest™ HBsAg

HBsAg Lateral Flow

The AllTest™ HBsAg Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Hepatitis B Surface Antigen (HBsAg) Virus in serum, plasma or whole blood. For professional in vitro diagnostic use only.

PERFORMANCE CHARACTERISTICS:

- ✓ Direct Testing from serum, plasma or capillary whole blood.
- 🕒 Reliable Results in 15 minutes.
- ✓ Sample: 3 drops of serum or plasma / 3 drops of wb or fingerstick.
- ✓ Sensitivity: 99.9%
- ✓ Specificity: 99.9%

Optional: Sterile lancet and Alcohol pad.

