

HCG Rapid Test HIGH SENSITIVITY

The HCG H.S. detects the presence of hCG in urine and serum specimens at the sensitivity of 10 mIU/mL and is an aid in the early detection of pregnancy.



- Simple, Accurate, Fast.
- Direct Testing on serum and urine.
 - ⌚ Reliable Results in 5 minutes.
- Sensitive: detects as low as 10 mIU/ml HCG
- Sensitivity : >99.9% code: 152000
- Specificity: >99.9% code: 152000
- Sensitivity : 100.0% code: FHC-U202
- Specificity: 100.0% code: FHC-U202

ORDERING INFORMATION:

Product Code
152000
FHC-U202

Product
HCG Rapid Test
HCG Rapid Test

Pack-Size
30 test
25 test

HCG Rapid Test ULTRA SENSIBILE

Test rapido immunocromatografico Lateral Flow ad alta sensibilità per la rilevazione qualitativa della gonadotropina corionica umana (hCG) in campioni di urina o siero.

- Semplice, Accurato, Veloce.
- Test diretto su siero e urina.
- ⌚ Risultato in 5 minuti.
- Limite di Rilevazione: 10 mIU/ml
- Sensibilità : >99.9% codice 152000
- Specificità : >99.9% codice 152000
- Sensibilità : 100.0% codice: FHC-U202
- Specificità: 100.0% codice: FHC-U202

