

HEPATITIS DISEASE

HCV

IMMUNOASSAY LATERAL FLOW

AllTest™ HCV

HCV Lateral Flow

Test rapido qualitativo di screening per la determinazione del virus dell'Epatite C Anticorpi (HCV) in campioni di siero, plasma e sangue intero. Test ad esclusivo uso professionale in vitro.

SPECIFICHE:

- ✓ Test diretto su siero e plasma e sangue intero e capillare.
- 🕒 Tempo di esecuzione: 10 min.
- ✓ Campione: 25 µL. siero o plasma / 50 µL. di sangue intero o da polpastrello.
- ✓ Sensibilità: 100%
- ✓ Specificità: 100%


Lancette sterili e Alcohol pad possono essere forniti separatamente.

INFO PER ORDINI:

REF
IHC-402

Product
HCV Rapid Test

Pack-size
25 test



AllTest™ HCV

HCV Lateral Flow

The AllTest™ HCV Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of antibody to Hepatitis C Virus in serum, plasma or whole blood. For professional in vitro diagnostic use only.

PERFORMANCE CHARACTERISTICS:

- ✓ Direct Testing from serum, plasma or capillary whole blood.
- 🕒 Reliable Results in 10 minutes.
- ✓ Sample: 25 µL. of serum or plasma / 50 µL. wb or fingerstick.
- ✓ Sensitivity: 100%
- ✓ Specificity: 100%

Optional: Sterile lancet and Alcohol pad.

CE
2934