

# HEPATITIS DISEASE

## IMMUNOASSAY LATERAL FLOW

**Biopanda**  
Reagents

## BIOPANDA® EPATITE E

**Nuovo Test Rapido** Lateral Flow per lo screening rapido qualitativo degli anticorpi del Virus dell'Epatite E in campioni di siero o plasma. Test professionale per esclusivo uso diagnostico in vitro.

## HEV IgG/IgM

### SPECIFICHE:

- Test diretti su siero o plasma.
- Differenziazione di IgG e IgM
- 🕒 Tempo di esecuzione: 10 min.
- Campione: 25 µL.
- Confezione: 20 test
- Sensibilità relativa: > 99,9%
- Specificità relativa: 100%
- Accuratezza analitica: 99,6%



## BIOPANDA® HEPATITIS E

The Biopanda HEV Rapid Test detects and differentiates between IgG and IgM antibodies to Hepatitis E virus (HEV) in human serum or plasma samples.

## HEV IgG/IgM

### PERFORMANCE CHARACTERISTICS:

- Direct Testing from serum or plasma.
- Detects and differentiates between IgG and IgM
- 🕒 Reliable Results in 10 minutes.
- Sample: 25 µL.
- Packaging: 20 test
- Relative Sensitivity: > 99,9%
- Relative Specificity: 100%
- Accuracy: 99,6%

### ORDERING INFO:

**REF. RAPG-HEV-001**  
**Product: HEV IgG/IgM**  
**Packaging: 20 test**

