

HEPATITIS DISEASE

IMMUNOASSAY LATERAL FLOW

BIOPANDA® HEPATITIS B



New Biopanda Hepatitis Rapid Tests are based on lateral flow immuno-chromatography to detect the presence of antigens or antibodies in human serum or plasma.

HBsAb, HBeAb, HBeAg, HBcAb
R.U.O.

COMMON CHARACTERISTICS:

- Direct Testing from serum or plasma.
- 🕒 Reliable Results in 15 minutes.
- Sample: 75 µL.
- Packaging: 20 test

PERFORMANCE CHARACTERISTICS:

HBsAb

- Relative Sensitivity: >99,9%
- Relative Specificity: 99,4%
- Accuracy: 99,6%

HBeAb

- Relative Sensitivity: 97,3%
- Relative Specificity: 97,9%
- Accuracy: 97,7%

HBeAg

- Relative Sensitivity: 99,9%
- Relative Specificity: 98,8%
- Accuracy: 99,1%

HBcAb

- Relative Sensitivity: 97,8%
- Relative Specificity: 97,7%
- Accuracy: 97,8%

ORDERING INFORMATION:

REF. RAPG-HBsAb-001 HBsAb 20 test

REF. RAPG-HBeAb-001 HBeAb 20 test

REF. RAPG-HBeAg-001 HBeAg 20 test

REF. RAPG-HBcAb-001 HBcAb 20 test

R.U.O.