




# Medizym<sup>®</sup> TRA human

**New Enzyme immunoassay for the determination of autoantibodies to TSH receptor**



# Medizym<sup>®</sup> TRA human

- ▶ human recombinant receptor preparation
- ▶ human monoclonal detection antibody
- ▶ robust assay
  - ◆ stable at incubation temperatures up to 30°C
  - ◆ low influence of shaking speed
  - ◆ low influence of water quality
- ▶ Comparison of 133 serum samples (100 TRA positive patients, 33 healthy donors) resulted in 100 % agreement regarding positive/negative discrimination



- ▶ grey zone 1.0 - 1.5 IU/l  
positive > 1.5 IU/l

## Order information

Medizym<sup>®</sup> TRA human Order no. 3505

Meridian Healthcare<sup>®</sup>

MERIDIAN HEALTHCARE srl

Catania, ITALY

PHONE +39 095 725 68 69 | Fax: +39 095 725 4454

e-mail: [info@meridianhealthcare.it](mailto:info@meridianhealthcare.it)

[www.meridianhealthcare.it](http://www.meridianhealthcare.it)

