

MYOGLOBIN Rapid Test

The Meridian Healthcare Myoglobin Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of Myoglobin in human whole blood, serum, and plasma samples. This test is a tool to assist in the diagnosis of myocardial infarction (MI).

- Simple, Accurate, Fast.
- Sample: serum, plasma and W.B.
- 🕒 Reliable Results in 10 minutes.
- Sensitivity : 99.9%
- Specificity: 97.2%
- Cut-Off: 50 ng/ml.



MIOGLOBINA Test Rapido

Test rapido immunocromatografico a flusso laterale per la determinazione della Mioglobina nel siero, plasma o sangue intero umano.

- Semplice, Accurato, Veloce.
- Test su siero, plasma e sangue intero.
- 🕒 Risultato in 10 minuti.
- Sensibilità : 99.9%
- Specificità : 97.2%
- Cut-Off: 50 ng/ml.



EUROPE Ordering Information:
Myoglobin RAPID TEST
Product Code: **CMY-402**
Packaging: **25 test**

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)