


## NADAL™ CELIAC DISEASE tTg IgA Rapid Test

The **Celiac Disease tTG** rapid test is an immunochromatographic test designed to detect, in human blood, IgA-type antibodies against human tissue transglutaminase (tTG), the main autoantigen recognized by the anti-endomysial antibodies.



- Direct Testing in whole blood.
- Sample: one drop of blood ( 10µL. )
- 🕒 Reliable Results in 10 minutes.
  - 10 Sterile lancets
  - Sensitivity tTG : 95,8%
  - Specificity tTG: 96,6%
- Long shelf life of 28 months



## NADAL™ CELIAC DISEASE tTg IgA Test Rapido

Il dispositivo **Celiachia NADAL™ tTG** è un test rapido Lateral Flow per la determinazione degli anticorpi anti-Tranglutaminasi IgA su campioni di sangue intero. L'utilizzo è particolarmente indicato nel caso di pazienti pediatrici.

- Semplice, Accurato, Veloce.
- Test diretto su sangue intero.
- Campione: 10µL. di sangue venoso o capillare.
- 🕒 Risultato in 10 minuti.
- Sensibilità tTG : 95,8%
- Specificità tTG: 96,6%
- Kit completo di lancette pungidito.



### ORDERING INFORMATION

**Codice Dispositivo**  
790004N-10

**Dispositivo**  
NADAL™ Celiac tTG IgA Rapid Test

**Confezione**  
10 test