

N-Terminal Pro-B-Type Natriuretic Peptide

NT-proBNP Lateral Flow test applies colloidal gold immuno-chromatography to detect N-terminal B-type natriuretic peptide precursor (NT-proBNP) in serum or plasma samples semi-quantitatively with standard colorimetric card. This test is used as an aid in the clinical diagnosis, prognosis and evaluation of Heart Failure (HF).

FEATURES:

- Simple, Accurate, Fast.
- Sample: serum or plasma.
- 🕒 Reliable Results in 15 minutes.
- Limit of detection: 200 pg/mL.
- Semiquantitative.



N-Terminal Pro-B-Type Natriuretic Peptide

NT-proBNP è un test rapido immunocromatografico a flusso laterale per la determinazione semiquantitativa dell'NT-proBNP su campioni di siero o plasma. Il test è di ausilio per diagnosticare e valutare la gravità dello scompenso cardiaco.

SPECIFICHE:

- Saggio semiquantitativo.
- Campione: siero o plasma.
- 🕒 Risultato in 15 minuti.
- Limite di rilevazione: 200 pg/mL.
- Kit provvisto di scala colorimetrica.

I risultati sono riportati in cinque intervalli di concentrazione:

- minore di 200 pg/ml
- tra 200 e 450 pg/ml
- tra 450 e 900 pg/ml
- tra 900 e 4000 pg/ml
- maggiore di 4000 pg/ml.



EUROPE Ordering Information:

NT-proBNP RAPID TEST

Product Code: HGCG03

Packaging: 20 test

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)