

Anti-Mullerian Hormone Elisa Test

This Test kit is used for in-vitro quantitative determination of the concentration of anti-Mullerian hormone in human serum sample, mainly used for clinical diagnosis of ovarian reserve function.




- Quantitative ELISA test
 - Sample: Serum
 - 🕒 Incubation Time: 80 min.
 - Volume Serum: 50 µl
- Reference Range Serum: 0.1 ng/ml - 17.5 ng/ml
 - Controls: Two Level Dried
 - 96 test/kit
 - CE marked

Anti-Mullerian Hormone Elisa Test

Test ELISA competitivo in fase solida per la rilevazione quantitativa dell'ormone AMH in campioni di siero.

- Test ELISA enzimatico competitivo.
- Tempo del saggio: 80 min.
- Volume campione Siero: 50 µl
- Linearità: 0.1 ng/ml - 17.5 ng/ml
- Controlli: 2 livelli inclusi
- Numero test: 96
- CE marked



EUROPE Ordering Information:

AMH AccuBind ELISA

Product Code: 9725-300A

Packaging: 96 test



Meridian Healthcare®

Meridian Healthcare srl

68 - Tremestieri Etneo

Catania, IT 95030

Phone: 095 - 725 68 69 (0039)