

PTH Intact ELISA KIT



Parathyroid Hormone (PTH) ELISA
Intact-PTH: For the Diagnosis of Primary
Hyperparathyroidism and Management of Chronic Renal
Failure & Renal Bone Disease .



Monobind Inc.
Lake Forest, CA 92630, USA

AccuBind
ELISA Microwells

elisa PTH intact

Parathyroid Hormone (PTH) ELISA
determinazione quantitativa del PTH intatto (ormone
paratiroideo) nel siero umano.

PTH Intatto

ELISA KIT Paratormone



SPECIFICHE

- ❖ Fase solida: enzyme-linked immunosorbent assay - ELISA
- ❖ Metodo: biotina-streptavidina
- ❖ Reazione: colorimetrica direttamente proporzionale
- ❖ Controlli: 1 alto e 1 basso inclusi nel kit
- ❖ Campione: siero
- ❖ Sensibilità: 0,49 pg/mL.
- ❖ Accuratezza: 95%
- ❖ Tempo di esecuzione saggio: 1 ora e 5 minuti
- ❖ CE Marked

Intact-PTH

ELISA KIT Parathyroid Hormone

FEATURES

- ❖ Name: PTH ELISA
- ❖ Description: ELISA for the quantitative determination of PTH Intact in human serum
- ❖ Design: ELISA
- ❖ Sample Volume: 25 µl serum
- ❖ Total Assay Time: 105 minutes
- ❖ Sensitivity: 0,49 pg/mL.
- ❖ Storage: 2 - 8 °C
- ❖ Kit: 96 Determinations
- ❖ CE Mark: CE Marked

Intact-PTH

ORDERING INFORMATION:

PTH INTATTO ELISA

REF. 9025-300

96 TEST



related products:

25-OH-Vitamin D (total) ELISA



CALCITONINA U.S. ELISA



CALCITONINA FAST ELISA

