

ROTA/ADENO COMBO & single Rapid Test

The StrongStep® Rota/Adeno Ag Rapid Test is a lateral flow immunoassay for the qualitative detection and differentiation of rotavirus and adenovirus antigens in fecal specimens. AVAILABLE: **ROTA/ADENO COMBO** - **ROTAVIRUS** and **ADENOVIRUS** single test.

FEATURES

- ⌚ Simple, Accurate, Fast.
- Sample: human feces.
- Ambulatory can be used.
- ⌚ Reliable Results in 15 minutes.
- Sensitivity Rota: >99,0%
- Specificity Rota: 99,5%
- Sensitivity Adeno: >99,0%
- Specificity Adeno: 99,4%
- CE marked



ROTA/ADENO COMBO & test singolo Rapid Test

Rotavirus e adenovirus sono i principali responsabili delle gastroenteriti. Rota-Adeno è un test rapido in grado di rilevare antigeni specifici nelle feci umane. Fornisce una risposta accurata, di primo livello, ed è di ausilio per il contenimento degli outbreaks, per evitare un trattamento antibiotico non necessario e per ridurre la mortalità infantile nei paesi in via di sviluppo. DISPONIBILI: **ROTA/ADENO COMBO** - **ROTAVIRUS** e **ADENOVIRUS** test singoli.

SPECIFICHE

- Semplice, Accurato, Veloce.
- Test diretto su campioni di feci.
- Utilizzo previsto anche in ambulatorio.
- ⌚ Risultato in 15 minuti.
- Sensibilità Rota: >99,0%
- Specificità Rota: 99,5%
- Sensibilità Adeno: >99,0%
- Specificità Adeno: 99,4%
- CE marked



EUROPE Ordering Information:



PRODUCT

ROTA+ADENO COMBO Ag
ROTAVIRUS Ag
ADENOVIRUS Ag

REF.

R0196C
IRO-602
IAD-602

N. Test

25 card
25 card
25 card