

# RESPIRATORY SYNCYTIAL VIRUS Rapid Test

RSV Card test is a one step coloured chromatographic immunoassay for the qualitative detection of Respiratory Syncytial Virus antigens. It can be used directly with nasal swabs or nasal wash or nasal aspirate specimens.


## FEATURES

- Simple, Accurate, Fast.
- Nasal Wash or Aspirate specimen.
- 🕒 Reliable Results in 15 minutes.
- Sensitivity: 92.7%
- Specificity: 98.5%
- Accuracy: 95.6%

### ORDERING INFORMATION:



Code	Product	Pack-size
IRS-502	RSV Rapid test	20 test



# VIRUS RESPIRATORIO SINCIZIALE Rapid Test

Test rapido lateral flow per la rilevazione qualitativa del Virus Respiratorio Sinciziale in secrezioni nasofaringee e aspirato nasale.

## SPECIFICHE

- Semplice, Accurato, Veloce.
- Tampone o Aspirazione Nasofaringea.
- 🕒 Risultato in 15 minuti.
- Sensibilità: 92.7%
- Specificità: 98.5%
- Accuratezza: 95.6%
- Per Uso Diagnostico in Vitro

