

Zinc – (5-Br-PAPS) MONOREAGENT

Colorimetric test for the determination of Zinc in human Serum, Plasma and Urine. Monoreagent END-POINT. For in vitro diagnostic use only.



- Colorimetric End-Point Assay.
- Single Reagent Liquid Stable Test Kit .
 - 🕒 Reliable Results in 5 minutes.
- **Linearity: Up to 500µg/dl (61,2µmol/l)**
- Sample: Serum, Plasma and URINE .
 - Standard included.
- Applications are available for many analyzers after request.



ORDERING INFORMATION:



Product Code
ZF01000050

Product Open System
Zinc - Zinco Monoreagente

Pack-Size
2x50 ml.

ZINCO – FLUID MONOREAGENTE

Saggio colorimetrico END-POINT per la determinazione dello Zinco su campioni di siero, plasma e URINE.

- Metodo Colorimetro END-POINT .
- Monoreagente Liquido pronto all'uso.
- Calibratore incluso.
- 🕒 Risultato in 5 minuti.
- Saggio su siero, plasma e URINE.
- Automatizzabile su tutte le piattaforme open system.
- **Linearità: 500µg/dl (61,2µmol/l)**

DEDICATED PACKAGING

Product Code
ZF01SEL025
ZF01911050

ANALYZER
SELECTRA Series
HITACHI 717
HITACHI 902
HITACHI 911

Pack-Size
5x25 ml.
4x50 ml.