

Cardiac Troponin I (cTnI) Rapid Test

The Cardiac Troponin I Rapid Test is a lateral flow chromatographic immunoassay for qualitative detection of Cardiac Troponin I (cTnI) in whole blood, serum, and plasma samples. This test is a tool to assist in the diagnosis of myocardial infarction.

- Simple, Accurate, Fast.
- Sample: serum, plasma and W.B.
- 🕒 Reliable Results in 10 minutes.
- Sensitivity : 99.4%
- Specificity: 99.0%
- Cut-Off: 0.5 ng/ml.



Troponin I (cTnI) Marker Test Rapido

Test rapido immunocromatografico qualitativo a flusso laterale per la determinazione della TROPONINA in campioni di siero, plasma e sangue intero.

- Semplice, Accurato, Veloce.
- Test su siero, plasma e sangue intero.
- 🕒 Risultato in 10 minuti.
- Sensibilità : 99.4%
- Specificità : 99.0%
- Cut-Off: 0.5 ng/ml.



EUROPE Ordering Information:
Troponin I (cTnI) RAPID TEST
Product Code: CTI-R402SQ
Packaging: 25 test

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)