

ENTAMOEBA + GIARDIA + CRYPTO Rapid Test

The *Entamoeba+Giardia+Crypto* Combo Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of *Entamoeba histolytica* antigens, *Giardia lamblia* and *Cryptosporidium* antigens in human feces specimen.

Product Features

- Simple, Accurate, Fast
- Direct Testing on Faeces
- Ambulatory can be used
- 🕒 Reliable Results in 10 minutes
- Entamoeba Sensitivity: 97,7%
- Entamoeba Specificity: 99,2%
- Giardia Sensitivity: 94,5%
- Giardia Specificity: 97,6%
- Crypto Sensitivity: 96,7%
- Crypto Specificity: 98,0%
- CE marked



ENTAMOEBA + GIARDIA + CRYPTO Test Rapido

Test rapido immunocromatografico Lateral Flow per la rilevazione simultanea e differenziata degli antigeni Entamoeba Histolytica/Giardia Lamblia/Cryptosporidium in campioni fecali.

Specifiche

- Semplice, Accurato, Veloce
- Test diretto su campioni fecali
- Utilizzo in laboratorio e in ambulatorio
- 🕒 Risultato in 10 minuti
- Entamoeba Sensibilità: 97,7%
- Entamoeba Specificità: 99,2%
- Giardia Sensibilità: 94,5%
- Giardia Specificità: 97,6%
- Crypto Sensibilità: 96,7%
- Crypto Specificità: 98,0%
- CE marked



EUROPE Ordering Information:
Entamoeba/Giardia/Crypto Combo Rapid Test
Product Code: IEGC-635
Packaging: 10 test

