


NORO+ROTA+ADENO+ASTROVIRUS Rapid Test

The Norovirus/Rotavirus/Adenovirus/Astrovirus Combo Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of norovirus, rotavirus, adenovirus and astrovirus in human feces specimens to aid in the diagnosis of norovirus, rotavirus or adenovirus infection.

- ❖ Simple, Accurate, Fast.
- ❖ Direct Testing on faeces
- ❖ Simultaneous detection of NORO+ROTA+ADENO+ASTRO
- 🕒 Reliable Results in 15 minutes.
- ❖ NORO Sensitivity : 95,7%
- ❖ NORO Specificity: 91,7%
- ❖ ROTA Sensitivity: 97,3%
- ❖ ROTA Specificity: 97,1%
- ❖ ADENO Sensitivity: 95,2%
- ❖ ADENO Specificity: 97,7%
- ❖ ASTRO Sensitivity: 97,0%
- ❖ ASTRO Specificity: 97,2%



NORO+ROTA+ADENO+ASTROVIRUS Test Rapido

Test rapido qualitativo a lettura visiva per la rilevazione differenziata di ADENOVIRUS/ROTAVIRUS/NOROVIRUS/ASTROVIRUS Antigeni su campioni fecali. Il NORO utilizza anticorpi monoclonali specifici di Genogruppo 1 e Genogruppo 2.

- ❖ Semplice, Accurato, Veloce.
- ❖ Test Diretto su campioni fecali.
- ❖ Rilevazione simultanea di NORO+ROTA+ADENO+ASTRO
- 🕒 Risultato in 15 minuti.
- ❖ NORO Sensibilità : 95,7%
- ❖ NORO Specificità: 91,7%
- ❖ ROTA Sensibilità: 97,3%
- ❖ ROTA Specificità: 97,1%
- ❖ ADENO Sensibilità: 95,2%
- ❖ ADENO Specificità: 97,7%
- ❖ ASTRO Sensibilità: 97,0%
- ❖ ASTRO Specificità: 97,2%



ORDERING INFORMATION:

Code No.
IMVD-645

Product
NORO+ROTA+ADENO+ASTRO Rapid Test

Packaging
10 Test