


# ELECTROPHORESIS AUTOMATED SYSTEM

## ADALYA IFE




**SELEO**  
SRL

Meridian  Healthcare®

**MADE IN  
ITALY**

# ELECTROPHORESIS AUTOMATED SYSTEM



## VERSATILE

Allows to perform the most common techniques: Serum proteins, Lipoproteins, Urinary Proteins, Immunofixation.

## COMPLETE AUTOMATION

AdalyaIFE automatically arranges the loading of all the reagents, to standardize the delicate electrophoretic process.



Serum Proteins



Haemoglobins



Lipoproteins



Integrated webcam for reporting of results



Automated Ag/Ab reaction contact

## PRODUCTIVITY

Serum-Proteins : 40 samples  
Lipoproteins: 40 samples  
Haemoglobins: 40 samples  
Immunofixation: 5 samples

## SOFTWARE FEATURES

User interface: monitor, keyboard and mouse

Minima recognition: automatic

Managing techniques: adjustable techniques in times, migration voltage/current

Operating system:



# ADALYA IFE



**SELEO** Meridian Healthcare®

DISTRIBUITO IN ITALIA DA

**MERIDIAN HEALTHCARE** srl

Via Caronda, 446 SC/A - 95100 CATANIA - ITALIA

Tel.: +39 095 7256869 Fax: +39 095 7254454

e-mail: [info@meridianhealthcare.it](mailto:info@meridianhealthcare.it)

web: [www.meridianhealthcare.it](http://www.meridianhealthcare.it)

## Your partner in Protein Electrophoresis