


ELECTROPHORESIS AUTOMATED SYSTEM

ADALYA 40




SELEO
SRL

Meridian Healthcare®

MADE IN
ITALY

ELECTROPHORESIS AUTOMATED SYSTEM



FAST/VERSATILE

4/40 samples, up to 5
electrophoresis membranes
8 samples each.


Allows to perform the most common
techniques:

- Serum Proteins
- Haemoglobins
- Lipoproteins
- Urinary Proteins

COMPLETE AUTOMATION

Adalya automatically arranges
the loading of all the reagents
to standardize the delicate
electrophoretic process.

During test all the liquid are
constantly measured and
adjusted, minimising the
variables introduced by reagents
level variations and
to ensure the test preparation in true
complete automation.



PRODUCTIVITY

Serum-Proteins : 24 samples/ 1h:15'

Lipoproteins: 16 samples/ 1h:15'

Haemoglobins: 16 samples/ 1h:15'


SOFTWARE FEATURES

User interface: monitor, keyboard and mouse

Minima recognition: automatic

Managing techniques: adjustable techniques in times, migration voltage/current

Operating system:



ADALYA 40



SELEO SRL Meridian Healthcare®

Your partner in Protein Electrophoresis

DISTRIBUITO IN ITALIA DA

MERIDIAN HEALTHCARE srl

Via Caronda, 446 SC/A - 95100 CATANIA - ITALIA

Tel.: +39 095 7256869 Fax: +39 095 7254454

e-mail: info@meridianhealthcare.it

web: www.meridianhealthcare.it